

Product datasheet for SC212526

FAM43B (NM_207334) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FAM43B (NM_207334) Human 3' UTR Clone
Symbol:	FAM43B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_207334
Insert Size:	1087 bp
Insert Sequence:	>SC212526 3'UTR clone of NM_207334 The sequence shown below is from the reference sequence of NM_207334. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CCCAGGGAGCGGGCGGGCCAGGCGCGCTGAGAGCCGAAGGACAGGACTCGCAGCCCCAGGCCCGACCCG
CCAGACTCACAGCCTCCAACCCCGGCCCTGCCCGCTTCGGCTGCCCGGCCCGGCCCGTGTCTCCCC
CGTGGTCTCCGTGTTGTCCGCCCGCCGCTCATTCTGGCTCAGGGTGATGCCTGATACGCCCTTGTT
ATTGGGGGTGTTCTCTCTCCCCACACCCGAGTTTCCCGGGCTGCCATTGTGGACCCGCCCTAT
GCTTTACCTAGTCTCTTTGCCACAGACCTCCTCATTCCCTCCAAAACATCCTCTCAAGAGAAGGG
AGGAGAAGTTTCAAGAAATCAGGAGGGGTGGTTTGGACCCTGGGCAGGGTGGAGGCAGTGACCTTGCC
CTTGGTCCCTCTAGCCTTCTTCCTGTGCAAAAAAATGACCCTGGAGAGGCATTCTTGTAAGGAGAAG
AATCTAGCGGCCGGGAGAATTGGGGCCGGCGCGGTGGGCAGAGTCCGCTGTATACACACAGGGA
GGAATTCTCAGCCCCAAGCCCCGCTCTCTACGCTTGGAGGACTCCTGTGACTTCACTGCTCTGCCTC
TGGAGAACTGGGAGAGTCCTACCGACGTTCAAACAACAGGTTAGGCCAGGTAACAGCCTGCACCAG
GCCGCTGCCACGCTCTGCCCTGGCACCCCCAGGGATTCTTGCCCATCCATCTCTCTGCAGACGG
ATGTGTGTGGCCCCCTCTAGGTGCCCCACAACAGGACCAAGATGGGGCTCCCAAAGGAGGTAAGGAG
AACCTTTGGCAGGTGCTTAGGACTGACTACCTAGAAAGTAGACGCAGCAGAGTTGCTCCCAAGTCGA
GGCTCCTCAGAGCAGGTGGGTCTGACAGCAGTGGATTCTCCAGCAGGATGAGGAAGGAGGGTGTGTT
AACCAACCAAGGGAGTGGGCCCCCACCAGGTGTCTCCGAAGACCACAAAAGCCCAAAGATCTATG
TGCTACTGATCATTGTAAATAAAGTGGACCTGCTTTACAGCCCTGTCTACTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_207334.3</u>
Locus ID:	163933
MW:	38.5