

Product datasheet for **SC212529**

C19orf42 (SMIM7) (NM_024104) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: C19orf42 (SMIM7) (NM_024104) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: SMIM7
Synonyms: C19orf42
ACCN: NM_024104
Insert Size: 1128 bp
Insert Sequence: >SC212529 3'UTR clone of NM_024104
 The sequence shown below is from the reference sequence of NM_024104. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCTGCATGATTGTGCTGTTTCGGCTCTTGAATCCAGCGATGAAACCAGGAACACTCTTCCCGGGATG
CCGAGTCTCCATTCCCTCATTCCCTGATGACTTCAAGAATGTTTTGACCAGAAAACCGACAACCTTCCC
AGAAAGTCCAAGCTCGTGGTGGGTGAAAAGTGTTCGCCGAGGTGTCATGGTTTTCCAGCCACGTCCC
TGTTTTCAAAGATAGTTTCACTTTGGTCTCTGAATTGAAATGCTGTCTACTGAAAGGGTTTCAGGAGCG
TTTATGTAAGGGGCTGTGATGAAATTGCATTCCCATAGATAAAAAGAAAATCATTTCTATCCAGAGAT
CTGAGCAGAAGGATTGGCTTGTAGTTAACACAGCCGATTTTTGGACATTCAGTGTACTTGCTGAG
TCTGACAGCCTCTGGGCCCGCCAGGGGCCCTGTTAACAACTGCTTTCACATCCCAACAGGGTCTGCT
TGGCCACTCAGTGCAGCTGCGATTAACCCTAAAGGCTTTAAGGAACGGGCCACCTGTAACAGAGACACC
AGCCTTCTGTATAGACTAAATTTAGCAAGAGTGTGAGCTAGTTCCTGGTGAAGTGTTCACACA
GAAGACATGTGGAGCAGTTGTGGGATATTAAGGGAACTTTCTCTGCCTTGACCCCTTTGTTAAATA
AAATGACTTTGGGAGCCATTCATTGTACAGTTGCAGGAATGAGAGTGATTTATGATGTGGTACATTGG
GACCATGTTCTAAAACCTTGGGTTTCTGAGTCTGCTTTTTGAGTAGGTGATTTTGAGGTTGAAAAACA
GGGCCTTCATCTAGGAAATACAGCATTTCAGAAAGCTTTTGAAGGGAACTCCTGGTTTTGTTGCC
AAAATGAAACGCCCGGGTTGGCGCTGAATCCACAACCTGTGTGATTTGCTTGTGAGTTTTTTGTTGT
CTGTTTTTTTTGTTTGTGTTTATACCAATAAGAATGAGCCTGAATGTTGGTGGTTTTTGAATCCTG
ACTTGGAGGTAACCTGGAGGAAGGAAAAAAGTAAATATGCAGGCTTTTAGGACTGAGTAGCCTTGA
AATAAACCTCATTCTAAAAAGGA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024104.4</u>
Locus ID:	79086
MW:	42.4