

Product datasheet for SC212550

ZNF563 (NM_145276) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF563 (NM_145276) Human 3' UTR Clone
Symbol:	ZNF563
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145276
Insert Size:	1147 bp
Insert Sequence:	>SC212550 3'UTR clone of NM_145276 The sequence shown below is from the reference sequence of NM_145276. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAAAGACTCACTGGAGAGAAACAATATGAATGTAACTATTATAAAGCCTTCTATATTTCCAGTTCCC
TTTGATATCATGAATAAACTTACACTGGACAGAAACCTATAAATGTAATGTGGGAAAGCCTTTCATA
ATTTTTATTCCTTTCAAAATCATGAAATATCTCACACTGGAGAGATGCTGTATGATATAAGGAATGTGG
GAAAGTGATCGATTTTTTTTTTAACACCATTTGTCAACATAAAACAACACAGAGCAGAAAACTCTGTG
AGTGTAATAATATGTAGGGAAGCCTTCAGTCATTTTGGTAACCTAAAAGTACATGAGGCTGGGCACAGTG
GCTCACGCCTGTAATCCCAGCACTTAGGGTGGGCAAATTGCTTGACCTCAGGAATTTGAGGCCTGCCTG
GGCAAAATGGTGAAACCATGTCTCTACAAAAAATAGAAAAAGTAGCAGGTGTGGTGGTATGCCCCTGTA
GTCCCAGCTACTTGGGAGGCTGAGGTGGGAAGATCTCTTGAGCCCAAGGGGTGGGGTTGCAGTGAGCC
AAGATGCACTGCACCCAGCCTGGGCGACAGAGCCAGACCCTGTCTCAAAAAAAAAAAAAAGTACATG
AAAGGATTCACCTCTGGAGAGAAGTCATATGAATGTAAGAAATGTAGAAAAGCATTGAGTTCTCTCAGGT
CCTTGATAGACATAAAAGGACTCACTGGAGAGATACTCTATAAATGTGTGAAATGTAGGGAAGTATTC
TTTAATGTTATTTTCAATTCAGATAGAAATAATTCTCATTGGAGATAAACCTATGAATGTAAACACATA
GTAAAGCCTTAATTTTAAGCCTCAGTTTCTTTCAAATACAGTTATCCCCTAATATTTGCAGGAGCGTGG
TTCTACCACCCCCTGGGCATACATAAATCGACACATGCCAAGTCCTTTTTTAAAAATGACATATTTGGT
ACAAAATTACCCGGATGTGATGGCAGATGCCTGTAATCCCAGCTACTCAGGAGGCTGAGGCAGGATAAT
CACTTGAATCCGGGAGGCGGAGGTTGCGGTGAGCCAAGATCGCGCCACTGCACCTCCAGCCTGGGCAGCA
AGAGCAAAACAACCTCTGTACACACACAAAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_145276.3</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	147837
MW:	43.6