

Product datasheet for SC212587

Kindlin 2 (FERMT2) (NM_001134999) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Kindlin 2 (FERMT2) (NM_001134999) Human 3' UTR Clone
Symbol:	Kindlin 2
Synonyms:	KIND2; mig-2; MIG2; PLEKHC1; UNC112; UNC112B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001134999
Insert Size:	1134 bp

Insert Sequence: >SC212587 3'UTR clone of NM_001134999
The sequence shown below is from the reference sequence of NM_001134999. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCTACAACTTACCAGTGGTTGGGTGTGAATAGGAATACTGTTTAAATGAACTCCACGGCCATAACAA
TATTTAACTTTAAAAGCTGTTTGTATATGCTGCTTAATAAAGTAAGCTTGAAATTTATCATTTTATCA
TGAAAACCTCTTTGCCTTACCAGACCAGTTAATATGTGCACTAAACAAGCAGACTATTAATCTATCAT
GTTATGATATAATAAACTTGAATTTGTACACATTCCTTAGGGCCATGAATTGAAAACCTGAAATAGTGG
GCAAAATCAGGAACAAACCATCACTGATTTACTGATTTAAGCTAGCCAAACTGTAAGAAACAAGCCATCT
ATTTTAAAGCTATCCAGGGCTTAACTATATGAACTCTATTTATCATGTCTAATGCATGTGATTTAATG
TATGTTTAAATTTGATATCATGTTTAAAATATCCTACTTCTGGTAGCCATTTAATTCCTCCCCCTACCC
CCAAATAAATCAGGCATGCAGGAGGCCTGATATTTAGTAATGTCATTGTGTTTGACCTTGAAGGAAAAT
GCTATTAGTCCGTCGTGCTTGATTTGTTTTGTCCTTGAATAAGCATGTTATGTATATTGTCTCGTGT
TTTATTTTACACCATATTGTATTACACTTTTAGTATTCACCAGCATAATCACTGTCTGCCTAAAATAT
GCAACTCTTTGCATTACAATATGAAGTAAAGTTCTATGAAGTATGCATTTTGTGTAACCTAATGTAATA
CACAAATTTTATAAAATTTGACAGTTTTTTAAAACCTACTCACAACCTAGCAGATGGCTTAAATGTAGCA
ATCTCTGCGTTAATTAATGCCTTTAAGAGATATAATTAACGTGCAGTTTTAATATCTACTAAATTAAG
AATGACTTCATTATGATCATGATTTGCCACAATGTCCTTAACTCTAATGCCTGGACTGGCCATGTTCTA
GTCTGTTGCGCTGTTACAATCTGTATTGGTGCTAGTCAGAAAATCCTAGCTCACATAGCCAAAAGGG
TGCGAGGGAGAGGTGGATTACCAGTATTGTTCAATAATCCATGGTTCAAAGACTGTATAATGCATTTT
ATTTTAAATAAAAAGCAAACTTTTATTTAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATCCACCGCCGCTTCTATGAAAGG
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001134999.2</u>
Summary:	Scaffolding protein that enhances integrin activation mediated by TLN1 and/or TLN2, but activates integrins only weakly by itself. Binds to membranes enriched in phosphoinositides. Enhances integrin-mediated cell adhesion onto the extracellular matrix and cell spreading; this requires both its ability to interact with integrins and with phospholipid membranes. Required for the assembly of focal adhesions. Participates in the connection between extracellular matrix adhesion sites and the actin cytoskeleton and also in the orchestration of actin assembly and cell shape modulation. Recruits FBLIM1 to focal adhesions. Plays a role in the TGFB1 and integrin signaling pathways. Stabilizes active CTNNB1 and plays a role in the regulation of transcription mediated by CTNNB1 and TCF7L2/TCF4 and in Wnt signaling. [UniProtKB/Swiss-Prot Function]
Locus ID:	10979
MW:	42.9