

Product datasheet for SC212614

NCCRP1 (NM_001001414) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NCCRP1 (NM_001001414) Human 3' UTR Clone
Symbol:	NCCRP1
Synonyms:	FBXO50; NCCRP-1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001001414
Insert Size:	1158 bp
Insert Sequence:	>SC212614 3'UTR clone of NM_001001414 The sequence shown below is from the reference sequence of NM_001001414. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGC
TCCTCCGTGTCTGTGCAGCTCCGGGAGTGACTGGCTGGCTCCTCTGTCCTGACCCACAGCACCTCCCT
GACCTTTAGGAGCCCCAACTCTTAGTCACCTCCTAGGCCTCTTATTTCTCCCTGGCCCTTGGCTTCTCA
CTTGATGGACAGCTTCACACACCCTTAAGCGGTGGACTCCAGCATTTTCCAGCACTGTCTGAGCCCCA
TGAGGGCGGAGCCACTCCTTGAAATTCAGTGCCCGACAGATGCTCTGGCACAGATGCTTTGTAGATCT
CTGTTGAGAGAATGCATAGACACCTGTGCCAAGGATGCTGAGGGCTGGTCTCTGCTTCTTTGAACCTC
ACTGAAACTGAATGCTCACTGCTGTGTTGCCAGCACACCAGCCAGGGCTGTGAACGGAGTGGGTGG
CAGCAATGTGTGTTGAAAGGGGAATGAAGCCATTCACTTCACTCAGTTCCTGTCCCATTTAACC GCCC
CGATCCTTGATCTTCCATTACCTTCACATCCCGGGGTCTTCTGAACTGACCTTGACCTCTGATCTCTT
CACACATCTCCCCTTAGCATCTCCACTTACCTACTTTTTTTTTTTTTTTTGTAGATGGCATCTCACTCTG
TCACCCAGGCTGGAGTGTAGTGACACGATCTCGACTCACTGCAATTTCCACCTCTCAGGTTCAAGTGGT
TCTCCTGCCTCAGCCTCTCAAGTAGCTGGGATTACAGGTGCACAGCACCTCCCCGACTAATTTTATA
TTTTTAGTAGAGACGGGATTTCCGATGTTGGCCAGGCTGGTCTCAAACCTCCTGACCTCAAGTGATCTG
CCACCTTGGCTTCCAAAGTGCTGGGATTGCAGACGTGAGTCACTGCGCCAGCCATTCCATGTCTCT
TAAGTCTCAGAATCTCCCTAGCTCCCTCCAGGTGTCTGCAGTGTTGTCCCCTCAAAGCTGTCCACA
CCCTCCTCCGAGGACCCTTTGTGTATCTCCTCCAGCTACCGCAGAGCCACAAACCCAGGCATCTATCA
AAGTCCCTCATTATGAGGGTGGTGAAGACACAGACTGCGACCAAGACAGAAATATGAAATGTGAATG
ACAGCGTCCCCGTGTGTGAATGTGGGATTAAGCATTATCAACCTCTAA
ACGCGTAAGCGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001001414.2</u>
Summary:	Promotes cell proliferation.[UniProtKB/Swiss-Prot Function]
Locus ID:	342897
MW:	41.8