

Product datasheet for **SC212639**

ZNF187 (ZSCAN26) (NM_001023560) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF187
Synonyms:	SRE-ZBP; SREZBP; ZNF187
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001023560
Insert Size:	1138 bp

Insert Sequence: >SC212639 3'UTR clone of NM_001023560
The sequence shown below is from the reference sequence of NM_001023560. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AGATATCACCAACAAGACAACTGGCTGATGAGGTGTTCTCTCCTTGTAGAACATCAGAGAAGGCACA
TTGACTAGCAAACAGCACTTTAGGAAAAGTCACCGTAGCCCACTGTGGCATCAGAAAATCTTGGGGGC
TGAGTTGGAGGCTCCCTGCCTCTATTCTCTCCTTTGCTTTCTTGAAGTCAGCTTTGGACCACAATA
ATTTCACTGTAGATGATATGCTAGGATCAAAGTTAAACAGCATTCTTCACTGCAGGACATCTCAGAGCA
TGTAACATAACTGCATGATTATATACTCTAAGCAATAGAGAGCTTCATGACTGAGTAAGAGTTTTGAAG
TCAGCAGTGAATCAAGTGCCACAGATTTGCAGGCTTAAGCAGAACAAAGGAAGATTGATATTTTGAAG
TATGCTATAGCAGCTTTCTCTATGAAATAAACTGATGATGTTTGAAGTATACTACTCTCAAAGGTG
TCTTTAAAGTACAGGTTAATGGTGAACATTTCCCCAGTGGCTTCACCTCATTCTCCCACTGGCCTT
ACCCCTTCTTCCCACTGGAAGCATTTTCAAAGCAAAGATAATTTTCGCTGGTGAAGTTTCAGAACTT
GCCTTCAGGGTTAGCTTCATGTAATTTTACCATTTCCCATCCCCATTCCCACCACTTATGTCAAGAT
TCAAGTTATAAATTAACGTTTTACTTAGACTTTGAAAGAGATTTTCATGAGTAATTTGAATGAACCTTG
TGAATTAATCTGAATAGCAACTGTCTGACTTGATAACTCCAGTGCCAGTATAGTGGCTGCTGCATAGAC
ACTATTCAGTAAATGTCTGATGAAAGAAATGGACTTTTTCCCTTTAATACTGATGATGCAATTTTGAA
ATGTTGCTTGTCTGGAATACTTCGGAAGTTAAACAATAAAGTCAGCTTGGAGAGGAGATCATGATCTT
TATACTGTGAATAGCAGAATGTCACAGATGAAGAATATTAATTAAGGGCGTCCACTCCAAACATAT
CACCAGAGTGACAACTTCAGCCCTCAGCCTCTGCAATCCATATATATCCTGGTCTTGGGAGTCCATAGA
ATACTGTTTCCTTCTAATAAAGGTTTCAAACAAA
ACGCGTAAGCGGCCCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```



Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001023560.4</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	7741
MW:	43.6