

Product datasheet for **SC212646**

DNAJC25 (NM_001015882) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DNAJC25 (NM_001015882) Human 3' UTR Clone
Symbol:	DNAJC25
Synonyms:	bA16L21.2.1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001015882
Insert Size:	1153 bp
Insert Sequence:	>SC212646 3'UTR clone of NM_001015882 The sequence shown below is from the reference sequence of NM_001015882. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGC
CCTGGGCGGTTAACATTGTGGATGACTGAAGATTGATGGAATGCTACTATGCCAAACCTTAATTGTGA
TATTATTTTCATAACTGAATTATTTAGAAATGTATCAATTGACTGCTGCTCAGCAGTAACATAAAATTC
CTCAAGTATTTGATTAACAGAATAATGTCAAAATTTAAACCTTCCCTTAAACCTTTATACATAAGACA
TTTATGATTGTTCAATTTTATAATCTATTTGTGGATTTTGTAAAAGATTTACATGAAGATTTATTA
GTTGCCATTTAAATTTTATATGTTTAGTTAAAGATTTGACATGAAAATTTATTAGATACCATTTAA
AAATTTTATGTGTATGTGTTTATTCTTTGAGAATGTTACCTTACTGTTTGTAAATAGTGCTACATTTT
CTGCTTTCATGGCCTCTGTATTTTCACAAAACACCAAAAACAGCATTTAATTATATTATCATGAGTGTG
TTTCTGGACACAACTTCTGCAGTAGAATGACCTAACATGTGCTTTTCATTGCAGTCATTATAGGATTG
AAATACGTTCAAATAACCTCTCTAGGAAAGTCTTCTGCTAGAATTTCTCCCTCTATTCAATTATAATA
TTCTTTGTTTTTAAAGCCAGTCAAATATAATAGTCTTAATAAGATCAGAACTCTCCAGGAGAGTGAGT
CTACCTCACGTCCTTGTAGGATGATCTTGATTATAGTCTTATTATAGGACTATAACTGTATTCTCAAC
ATTTCTCCAGAAAGGACCTTGTAAGGAGTCTTTGTACACAGTATTGGTTTTTCCCTCTCTCTTC
ACTTAAAAAAAATAGCAAGGCAGAAATAGTGTATTGAAAAGTTGTTTCATCTATTATGAAGTCCTTG
AGTGGTGAAAAATCCGTTGTACATGAGAACATTTCTATGCATTTAAGCCAGAAACGAGGTACATGGCTG
TGTGCTCTTGTGCAACCAATGAAATGTGTTTTACATGTGTGGCAGTGCAAGTAAATAACACATTATT
TGACTGAATCAGGCATGATACTGCACCAAAGTGTGGTACATATTCACGGTAGTAAATCAGTACCCCTG
TTAAAGGATTTATCCATTGCTTCATATTAATAAAATGTTTACAATATA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001015882.3</u>
Locus ID:	548645
MW:	44.9