

Product datasheet for **SC212652**

PLAGL1 (NM_002656) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PLAGL1 (NM_002656) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PLAGL1
Synonyms: LOT1; ZAC; ZAC1
ACCN: NM_002656
Insert Size: 1159 bp
Insert Sequence: >SC212652 3'UTR clone of NM_002656

The sequence shown below is from the reference sequence of NM_002656. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGCCTCATTTCCATCATGCATTCAGATATTGATTTTTAAAGTGTATTTTTTCGTATTCTGGAAGATGT
TTTAAGAAGCATTTTAAATGTCAGTTACAATATGAGAAAGATTTGAAAACGAGACTGGACTATGGCT
TATTCAGTGATGACTGGCTTGAGATGATAAGAGAATTCTCGAACTGCATGTATTGTGCCAATCTGTCTCT
GAGTGTTTCATGCTTTGTACCAAATTTAATGAACGCGTGTCTGTAATCAAAGTCAAATATTGTCATAA
CCAACATCCAAAATGACGGCTGCTATATAAAGTGTGTCATATGGAATTTAATCGTAAGCCATGATC
ATAATGTTAACTAAATAACTTTATGTGGCACTGCCTAGTAAGGGAAGTAAAGGTTGGATTCTC
CAAATCTGGGAGAAATTTCAAATAAAGAAAATAACCTTTATATGATATACTATGACTAGGCTGTGTATT
TCTTTTCAGGGATTTTCTACCTTCAGGGTTGGATGTAGTTTAGTTACTATTACCATAGCCAACCTGTA
GTTTTACATATACATTTTCTGTGGAGCAATAGAGTTCTCCATTTACAGAAGCATTTTAAATGTAGTT
TGAAATTTTTCCACAAGATGCTGCAATGTGAGTTATCACTTCATTTATCTTAAAGAAAGACTAAACTGG
TTGTCAGTTACATCTGACAGAAAAAAAAAAAAAAAAATCACTGTGAACCAGGTTAAGTGGTAAAATAATCC
AGGCGTCAGTCAAAGGCATTTTGTGCTGACTTTAATATTGATTATATTTTTAACAGGAATTTAAGAAAATA
TACTGGAATTAATAATATATATATTAACAAGAATTTCTTTGCTCTGTCTAGCTTAACTACTAC
TCAAGCTGCTTAAGTTCTTAAGTATTGTTGTAATCACCAATAAATAAGTGCATTTGTAATTCATCAGT
CATTATTAGCTTTTATTAAGAAGATTACGTTTTACAATGTAATAATCTCTTGAATTTGGTATCT
TATTAATGAGTTTTAAAGATGTAACCTAACCTTTTTAAAGTCCATTGTCTTATGTTTTAGAGGC
TTTTCCGTAACATATATCTTACATATAATAAACTTTTCAAATCTTGCAAAAAA
AGCGGACCGACTTACGCGTAAGCGCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCC
CAACCTGCCATCAGAGATTTTCGATTCCACCGCCGC
```

Restriction Sites: SgfI-RsrII



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_002656.3</u>
Summary:	This gene encodes a C2H2 zinc finger protein that functions as a suppressor of cell growth. This gene is often deleted or methylated and silenced in cancer cells. In addition, overexpression of this gene during fetal development is thought to be the causal factor for transient neonatal diabetes mellitus (TNDM). Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding two different protein isoforms. The P1 downstream promoter of this gene is imprinted, with preferential expression from the paternal allele in many tissues. [provided by RefSeq, Nov 2015]
Locus ID:	5325
MW:	45.6