

Product datasheet for **SC212670**

Estrogen Related Receptor beta (ESRRB) (NM_004452) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Estrogen Related Receptor beta (ESRRB) (NM_004452) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ESRRB
Synonyms: DFNB35; ERR2; ERRb; ERR beta-2; ERRbeta2; ESRL2; NR3B2
ACCN: NM_004452
Insert Size: 1146 bp
Insert Sequence: >SC212670 3'UTR clone of NM_004452
The sequence shown below is from the reference sequence of NM_004452. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGGATCGCC
CAGGGAAAGCATCTCTGGCTCACCATGTAACATCTGGCTTGAGCAAGTGGGTGTTCTGCACACCAGGC
AGCTGCACCTCACTGGATCTAGTGTGGCTGCGAGTGACCTCACTCAGAGCCCTCTAGCAGAGTGGGG
CGGAAGTCTGATGGTTGGTGTCCATGAGTGGAAGCTGCTTTTATACTTAAACTCAGATCACAACAG
GAAATGTGTCAAGTAAACAATGGAATCCATCCAATGGGAAAGTTCCTGGTACTGAAGGGTCCATTGGAC
ACTCAGAAAAGAAGTTCAGGGGCAACTTCTTAGCTGGAATCCTGGCCAGATGAGGACCCTCTCCGGGG
AAGGGAGAGGACTGACTTAGTGAAGGTGGTGAAGTGAGGAGAGTTAGGGGAACCTTCCCCAGTGGA
ACAGATCTCAAGTTTACCCTAACCTGCCATTTCTGAAAATCTGTAAAGAGGAAACAGCCTGTCTCAG
CTGTACTCTCATGATACAGGTCATTTGAAATGAACCAAGAAATAAAACATGAAAATCCAACCATGGAGA
AGGTGGTATGGCTGGGTTTTGTTGGTCCCCTTGCTTATACGTTCTAAAGTTTCCAGACTGGCTTTG
TCACTTTGTGAACCTCGTCATGTGTGAAAACCAATCTTTGCATATAGGGAACCTTCCCTGGGCCACACTTT
AAGAACCAAGTAAGAGGCTCTCAAGACTCCAGCAGAGTCGGGAGGCCATGGCAGCGCCTTAGAGGAGCT
GGAACCTGCACCCACCTGTGTGGTGGGGGGGCTCCTTTCCCATAGACTCTGCCCTCCCTCTGTGC
AGATGGAAGTGGCAGGGGAGGGTGACCAGCTTGTGACAAGAAGACTGAAGGGTCCAGAGTCCATGTCA
CGGAACAGCACCAAGAAAAGCACTATGTGAAAAGATTGTTTTATTTCTAATAATGATAATATGGCTG
GAATGGCTTCTTAAGATGTATATATTTTTTAAATGGCAGTTCGCCATTGCAGCATCACCTACTTGTAT
GTCTTTCTGCCTGTATATGTTCTCCAGAAACCCCATGTAATCAAATGCCCTAGGATGCTTCCAT
CCTGGTCCCATGTATCTGGAATCTAATAATAAGGAAAGGAA
ACGGGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_004452.4</u>
Summary:	This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development. [provided by RefSeq, Jul 2008]
Locus ID:	2103
MW:	44.2