

Product datasheet for SC212680

LRRC8E (NM_025061) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LRRC8E (NM_025061) Human 3' UTR Clone
Symbol:	LRRC8E
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_025061
Insert Size:	1133 bp
Insert Sequence:	>SC212680 3'UTR clone of NM_025061 The sequence shown below is from the reference sequence of NM_025061. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAAGTGCGGGACAAGATGGAGGAGGAATGAAGCTGGGGTGGGGCCGTTTTAGGTAGAGCCTTAAAAATG
CTTCTGCCCTGGAATCTCAACCATTGTCTTCCAAGATAGGAAGCCAAGTGGGTCCAGGCCAGGAGATGG
GGGGGGCGGGGACAGCTGTGTCATCTTTCTGGGGCCCAGGAGGATCTGGGCTGGTTGTCTGGGGAGAC
AGACAGGATGTTGTGGAGCTGGGGTGGAACTGGTATGGAGGGATTAACTCAGTCATGGCATTCTCCGA
CCAAAACCAACACCTGTGTCTCTGGCAGGCTGGCTGGCCTTGCTCCCATCCCTAGAACTGCTGCCTCTCC
CTGGATATTCCAGCTCAATTAGTGCCACATATGGGGGAAACGACACATCCCAAGTGGGATTTCACAACT
CCCCCTCCCCATGCAACAAAGCAACTTACTTCTGGAGTTCTCTCCAAGGAGAGGACACAGACACAGTT
GTTTGCTGTGTTATATGTTAGCTCCGAACAATGGTTCTCATTTGGCTAAGCATCAAAATCACCTAGGGA
GTCGGTGCAAAACAAATATCCAGTCCCCTCCCCTGAAACACTGACTCAGGAGGTTTGGTTGGGGGCC
AGGAGTCTGTTCTAAATATTCCAGGTAGTTCTGGTGCAGGTAAGTGGCCCTGAGACAGTATGTTGGGA
AATGCTGACGTAAAGGTATCAGGGCCGGGCGCTGTGGCTCATGACTATAATCCAGCTGTTTGAGAGGC
CAATGCAGGAGGATGGTTGAGCTCAGGAGTTCGAGATCAGCCTGGGTAACATAGCGAGACCCACCTCT
GCCAAAAAAAAAAAAATAGGGTATCCAAATATCTAGATTCTGATCCCTTTTGAGGTCTAGACCCCTT
GAGAACTGATGAAGCCAGGCACCTCCTTCCTCAGGAAAATGCTGGTGTACAAATACACACAAAGCTCT
TCAGGCAGCTGATAGATTTCCCCAGAGAGCTATTCAAGGACTTCCTAAGGTGGGTGGACTGCAGGGTT
AGGACACCTGCTATAGAGGTGACATTTTCCAAGGACAAGCAGGGACTTTGGTCTTGACTGTTCTCTGG
GTGCTTAATAATAATACTAATTATGCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_025061.6</u>
Summary:	This gene encodes a member of a small, conserved family of proteins with similar structure, including a string of extracellular leucine-rich repeats. A related protein was shown to be involved in B-cell development. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2012]
Locus ID:	80131
MW:	42.2