

Product datasheet for SC212685

PRELID3A (NM_001142405) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PRELID3A (NM_001142405) Human 3' UTR Clone
Symbol:	PRELID3A
Synonyms:	C18orf43; HFL-EDDG1; SLMO1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142405
Insert Size:	1148 bp
Insert Sequence:	>SC212685 3'UTR clone of NM_001142405 The sequence shown below is from the reference sequence of NM_001142405. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATTGAACACTCTGAAAGCGCTGTGAGCTAAGGAGGCTGTGCCTGTGCTTGTCATAAATGGTGATGCTT
CCCAGAAGAAAGGAAAGAAGGAGGTGAGGAGAGGACGAGCAAGGACGGGACGGAAAGCACGTCTGGCT
GCGAGGGCACACAGGGCCGGGAGGAGGATTAGGCTGGACTGATCGGGGTCGAGGTGACCCCGGGAC
CAGATGTCCCTCCCTCATGAGCGTGCAACAGGGAGCAGCCCCAGCCAGCCCCACCGCCGCGTCTGCG
TGGATGGAGGGCTCGAGGCCCTCGCCCGTGTCCCTTAGGGGCGGCGCCATCGTCTGCCTGCCGAGT
GGCCAGCCCTTCTCTCCCTGCAGGTGACCCGCGCCGCCCCGCGCCGCCCCAGCTGTCTCCAC
TTAGGCTCCTCCACAAACGCTGGCGGGTGTGCCGCCCTCCAGCCTCTCCGCCCCGTCCTCGGTCT
GCAGATACCAAGCAGCTCTCCTGGGCCGCGGCTGTCCGGGCTCACATCTGTGTGGCTCGGTG
GGCTCTGACTGCGGGGCTCACGGGGCTGGGGCCATGGCTTTCTCTCGGCCCCAAATGCCTTTCTG
GAGTGGCTAAGCCTCTCCTGACGAGGCTCCTCATGCGGGGAGAGGGTGGGGTCCAGGGTTTCCCGG
CCTTTCTGCAGGCTGTGCCAGCCCCCTCAGTCTGCCACTGCCGTGCCGGGGGACATCCTGCGCCAG
GCTGGAGGCCCTGCCTCGCTGCTTGCCAGCCTCACCTGCTGGGGTGGAGGCCATGGATCCCTCATAGT
ATCCAAGTTTGAGTTCTGTGTTTTTCGCCTTATAGTTTTGTGAGTGTGGAGGACAACGACCAAGGT
CTTTTCCCCACCATCTTTAAATCAAAGTTGTTTAAACAACATTAATTTAAAAAATACATAAAAGCAA
CCCATGCTTATGGTTAAAAAAGAGTACAGGAGTCTGGAGTCCGTTGAAGGGTAAAGGCCTCCCGAAGC
TGCCAGTCGGCTCTTTTCGAAAGCCTCCAGAAGTTTTCTGTGCTTATGAAAGCGTATGTGTGCCT
TTCTCTTCCCTTCTTTTATTAACACAGATAGGTTTATGCTGGA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001142405.2</u>
Summary:	In vitro, the TRIAP1:PRELID3A complex mediates the transfer of phosphatidic acid (PA) between liposomes and probably functions as a PA transporter across the mitochondrion intermembrane space. Phosphatidic acid import is required for cardiolipin (CL) synthesis in the mitochondrial inner membrane.[UniProtKB/Swiss-Prot Function]
Locus ID:	10650
MW:	41.2