

Product datasheet for **SC212693**

Amyloid Precursor Protein (APP) (NM_001136129) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Amyloid Precursor Protein (APP) (NM_001136129) Human 3' UTR Clone
Symbol:	Amyloid Precursor Protein
Synonyms:	AAA; ABETA; ABPP; AD1; alpha-sAPP; APPI; CTFgamma; CVAP; PN-II; PN2; preA4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001136129
Insert Size:	1150 bp
Insert Sequence:	>SC212693 3'UTR clone of NM_001136129 The sequence shown below is from the reference sequence of NM_001136129. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TACAAGTCTTTGAGCAGATGCAGAACTAGACCCCGCCACAGCAGCCTCTGAAGTTGGACAGCAAAAC
CATTGCTTCACTACCCATCGGTGCCATTTATAGAATAATGTGGGAAGAAACAAACCCGTTTTATGATT
TACTCATTATCGCCTTTGACAGCTGTGCTGTAAACACAAGTAGATGCCTGAACCTTGAATTAATCCACAC
ATCAGTAATGTATTCTATCTCTCTTTACATTTTGGTCTCTATACTACATTATTAATGGGTTTTGTGTAC
TGTAAGAATTTAGCTGTATCAAAGTAGTCATGAATAGATTCTCTCCTGATTATTTATCACATAGCCC
CTTAGCCAGTTGTATATTATTCTTGTGGTTTTGTGACCCAATTAAGTCTACTTTACATATGCTTTAAGA
ATCGATGGGGGATGCTTCATGTGAACGTGGGAGTTCAGCTGCTTCTTGCCTAAGTATTCCTTTCTCTG
ATCACTATGCATTTTAAAGTTAAACATTTTAAAGTATTTTCAAGATGCTTTAGAGAGATTTTTTTCCATG
ACTGCATTTTACTGTACAGATTGCTGCTTCTGCTATATTTGTGATATAGGAATTAAGAGGATACACACG
TTTGTTCCTTCGTGCCTGTTTTATGTGCACACATTAGGCATTGAGACTTCAAGCTTTTCTTTTTTTGTC
CACGTATCTTTGGGTCTTTGATAAAGAAAAGAATCCCTGTTCAATTGTAAGCACTTTTACGGGGCGGGTG
GGGAGGGGTGCTCTGCTGGTCTTCAATTACCAAGAATTTCCAAAACAATTTTCTGCAGGATGATTGTA
CAGAATCATTGCTTATGACATGATCGCTTCTACACTGTATTACATAAAATAAATAAAATAAATAAACC
CGGGCAAGACTTTTCTTTGAAGGATGACTACAGACATTAATAATCGAAGTAATTTTGGGTGGGAGAA
GAGGCAGATTTCAATTTTCTTAAACAGTCTGAAGTTTCATTTATGATACAAAAGAAGATGAAAATGGAA
GTGGCAATATAAGGGGATGAGGAAGGCATGCCTGGACAAACCTTCTTTTAAAGATGTGTCTTCAATTTG
TATAAAATGGTGTTTTTCATGTAATAAATACATTCTTGGAGGAGCA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```



[View online >](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001136129.3
Summary:	This gene encodes a cell surface receptor and transmembrane precursor protein that is cleaved by secretases to form a number of peptides. Some of these peptides are secreted and can bind to the acetyltransferase complex APBB1/TIP60 to promote transcriptional activation, while others form the protein basis of the amyloid plaques found in the brains of patients with Alzheimer disease. In addition, two of the peptides are antimicrobial peptides, having been shown to have bacteriocidal and antifungal activities. Mutations in this gene have been implicated in autosomal dominant Alzheimer disease and cerebroarterial amyloidosis (cerebral amyloid angiopathy). Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2014]
Locus ID:	351
MW:	44.5