

Product datasheet for **SC212703**

ZNF385 (ZNF385A) (NM_015481) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF385
Synonyms:	HZF; RZF; ZFP385; ZNF385
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_015481
Insert Size:	1150 bp

Insert Sequence: >SC212703 3'UTR clone of NM_015481
The sequence shown below is from the reference sequence of NM_015481. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CACGGACCCATCCTCTTCTCCCGTACTGACCTCAACCCTGAACCCCTCCCATTTCAACTCCCCACCTCC
AGCCGGGAGCCAGGCGTCCGGGCTCCAGCCCGCCCCCTCTCCCGGCACTCCCTAAATGATCTCTCCTT
CCCCCCCCCACCAGAGATACGGGGTTCCAGGAAAGGGGAGGGGTAGCGGGGAGGGGGGTTTCAGAA
GGGGGGGAACACCCAGATCTCAGGGAACCCCGCCCCCTGCCTTTCCCTCTCCCTAGAAAAGGGGGGG
CCGTCTCACCCCGAGCCCCCTTGAGACACCCCCCTCCCAAAAGCCATGTCCATCCAGCCCTTCCCC
CCAAACCTAGCACAAAACGGGGTTCACAAGCCATGGTCGGGGTCCGGGGGGACAGAAATGGATTTTCT
TGGCAATAAGCGGACTCTGGGACTCCGGCTCCCTACCCCAAAGTGAAGCGCTTCCGTGAACACCCCGT
CCTCCGTAGGGGAGGGGAGCAGGCGGATCCTGGGTCCCTCATAAGCACTTTGGTTTTACCGCCTGCA
ACCTCACTGTGCCCGCCCGCACCATGCCCTAGCCCCAGGTCTAGCCGGGCCATTGCAAGGGGGCAGCA
CTTGGGGGCATCTCCGGCACTTGGGTGGGACCAAGGAGATGCCACCATAGACCTTTCCCTCGCCTTCTT
CCTCCCTAGTCCGGGTTCCATTCTTTTACCAGCACCATCGCCCAAGGGGTACCGAGGGGGCAAGGG
GTGTCCAGTCCAAGCCACCCCGCCTCGCCTTCCGCAAAACTGTGAGCAAAAAGCAATAGAAGCCTCG
CCCCCGCCCTTCCCTTCCGAGGATTCGCGAGTCTGTAGCCTCCCGATCCAAGTTCTAGACCTCAT
GGCTGTCCCTCCACCACTCACTCCACTGCACAACTCGGCGGGGGTGTGACACCTCTCCCCACCC
CCGACTCCGTGGTTTCCGTATCGTCAACCTTCAGCCGCCGACCCGGGAGGGGTCTGGCCTACACTGGT
CTTCCCTTCCCATCAACTCTTTCTGCTTGACAATGTAGCAACCCAGGCCCCACCCACGGTCTCTCC
CTTTTTCTCTCCCTGACAATAAAGTCTGAATTTGTTCTGCCCTCC
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```



Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_015481.3</u>
Summary:	Zinc finger proteins, such as ZNF385A, are regulatory proteins that act as transcription factors, bind single- or double-stranded RNA, or interact with other proteins (Sharma et al., 2004 [PubMed 15527981]).[supplied by OMIM, Oct 2008]
Locus ID:	25946
MW:	40.8