

Product datasheet for SC212714

SLC39A13 (NM_001128225) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC39A13 (NM_001128225) Human 3' UTR Clone
Symbol:	SLC39A13
Synonyms:	EDSSPD3; LZT-Hs9; SCDEDS; ZIP13
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001128225
Insert Size:	1163 bp
Insert Sequence:	>SC212714 3'UTR clone of NM_001128225 The sequence shown below is from the reference sequence of NM_001128225. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATGGTGTCTGTTCTCGCTCTTCGTGGATTAACTTTCCCTGATGCCGACGCCCTGCCCTGCAGCAATA
AGATGCTCGGATTCACTCTGTGACCGCATATGTGAGAGGCAGAGAGGGCGAGTGGCTGCGAGAGAGAAT
GAGCCTCCCGCCAGACAGGAGGGAGGTGCGTGTGGATGTATGTGGTGTGCACATGTGGCCAGAGGTGTG
TGCGCGAGACCGACACTGTGATCCCTGTGCTGGGTCCGGGGCCAGTGTAGCGCCTGTCCCAGCCATG
CTGTGGTTACCTCTCCTTGCCGCCCTGTACCTTCACCTCCCGGAGTAAGCAGCGAGGAAGAGCAGCAC
TGGTCCCAAGCAGAGGCCTTGCCCTGCTGGGACCCCGGAGTGAGAGCAGCCCAAGGATCCAGGGTGCA
GGGAACTCCAGAGCTGCCCACTCCCACTGCCCTCAGCACACACAGTCCCAGGCGGCCTAGGGG
CCAAGGCTGGGGCGGCTTTGGTCCCTTTTCTGGCCCTTCCTTCCCCTTCTAAGCCAAAGAAAGGAG
AGGCAGGTGCTCCTGTACCCAGCCCCACTCAGCACTGACAGTCCCAGCTCCTAGTAGTGAGCTGGGA
GGCGCTTCCTAAGACCCTTTCTCAGGGCTGCCCTGGGAGCTCATTCTGGCCAAACAGCCCTGGCAGC
ACCAGCAGCTCTTGCCACCTCCAGCTGCCAAACAGCAGCCTGCCGGGCAGGGAGCAGCCCCAGGCCAGA
GAGGCCCTCCCGTCCAGCTCAGGGATGCTCCTGCCAGCAGGGGCCAGGACTCCTGGAGCAGGCACA
TAGTGAGCCCGGCAGCCCTGCCAGCTCAGGCCCTTTCTTCCCATTGAGGTTGGGGTAGGTGGGG
GCGGTGAGGGCTCCACGTTGTACGCGCTCAGGAATGTGCTCCGGCAGAGTGTGAAGCCATAATCCCCA
ACCATTTCCTTGCTGACGCCCAGGTAAGTGTGCTGAGGCTGAGGCTGAGGCTGAGGCTGAGGCTGAGGCT
CACCCTGGATGTTTTAGAAAGTGGCCATCGAGAGGTCTGGATGGTTTTATAGCACTTTGCTGTGATTC
CGTTTGTATCTGTAATATTTGTTCTATAGATAAGATAACAAATAAATATTATCCACATA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001128225.3</u>
Summary:	This gene encodes a member of the LIV-1 subfamily of the ZIP transporter family. The encoded transmembrane protein functions as a zinc transporter. Mutations in this gene have been associated with the spondylocheiro dysplastic form of Ehlers-Danlos syndrome. Alternate transcript variants have been found for this gene. [provided by RefSeq, Jan 2016]
Locus ID:	91252
MW:	41.5