

Product datasheet for SC212781

ILDR1 (NM_175924) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ILDR1 (NM_175924) Human 3' UTR Clone
Symbol:	ILDR1
Synonyms:	DFNB42; ILDR1alpha; ILDR1alpha'; ILDR1beta
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_175924
Insert Size:	1074 bp
Insert Sequence:	>SC212781 3'UTR clone of NM_175924 The sequence shown below is from the reference sequence of NM_175924. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TCTCATAGTGAAGGAGTGTGGTCATTAGTCACCAAGCACAGCACAACTTCTGTGGCTACTTCTCGGC
TCCTGTGTGTCATCAGCATCACCTAGGTTTCCAGCTGACTTGGGAAGTCAAGTCTGAGTCTAACAGTT
TTTGGCTTAGATTCTGAGAATCAAATAGAAGAATTTTAAATACAAGAGTTTGAGATTGGGTATAGTGGC
TCACACCTGTAATCACAGCACTTTGGGAGGCTAAGAATCACTTGAGACTAGGAGTTCAAGATCAGCCTG
GGAAACATAGTGAGACCCGCTCTACAAAAAATATAAAATAGTTAGGTGTGGTGGCATGCATTGT
AGTCCCAGCTACTCAAGAGGCTGAGGCAGAAGGATCCCTTGAGCCAGGAGTGCAAGGCTGCAATGAGC
CAGGATCCCATGATCACACATCTGTATTCCAGCCTGGGCAATAGAGCAAGTCCCCATTACTAAAAAAC
CAAAAGGCCAAAAACAAAAAGTTAGAGTTCGAGGAATTACCAACTGTAGTTTTAGCCTTGTTTCATG
CTCTCTTGCAATTTATATAATCTCTGACTTGTAATGGACCCTGACTGGAATGTGATCCCTCAGGAAGT
TAGTAGCCTGAGTCTTTCAGTAGACTACACTGCCAGAACCTGGCCATTCTCAAAATGAGAACTGGG
AATGTTTAAGAAGAAATCAACATGTTTCAGGAAAAGGAAATCTATGGAGTATTATAGGGACATTCCTCA
TGGGAATGTATCTTCTCCATGGCATGTCTTGAGGGTCTTTCTTGTAGGAGTTTATCTGCCAGCCC
ATAAATGGACTATTTATTGTAAGTGTAGAAAATCACAGAGAAGCAGTTTTGCACCAGCCTTATTCTGT
GCCTTGTTTTCTCTTGCTCTTTTTTTACCTGTATATCTAATTTATATTTTCATATATATTGTGATTG
ATTGAAGTCACTTTAAATCCTTCTTGGAATGACACAGTATATAAATAAGGAAGAAAGAAACATGCCA
AGCTGAGCATGCTGCTCCAAATAAATATCTGCTTTCCTA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_175924.4</u>
Summary:	This gene encodes a protein that contains an immunoglobulin-like domain. The encoded protein may function as a multimeric receptor at the cell surface. The expression of this gene may be a diagnostic marker for cancer progression. Alternatively spliced transcript variants encoding multiple protein isoforms have been observed for this gene. [provided by RefSeq, Dec 2010]
Locus ID:	286676
MW:	40.2