

Product datasheet for SC212782

ARRDC2 (NM_001025604) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ARRDC2 (NM_001025604) Human 3' UTR Clone
Symbol:	ARRDC2
Synonyms:	CLONE24945; PP2703
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001025604
Insert Size:	1166 bp
Insert Sequence:	>SC212782 3'UTR clone of NM_001025604 The sequence shown below is from the reference sequence of NM_001025604. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GACATGAGGCCGCGCTGCATGACTTGCTACGCGCACAGGGACCCCTCGAGGAACAAGGTTGCACACCA
GCTTTAGCCACCATGACTGTGGGAGTGGCTGGACCAAGGGCTGACCTCCCCGACTGCATCAAAGTTG
GGGAACCAAGTCTCAGAGTGAGGCGGGGCTTTTCGGATATCACATGGGACAGAGGAAGAGCCGGCTG
GAATCTGACTTACCTGGACCGCTGTCCTTGTGAGGCATTGAATGCCAGTGACAGTATCCGAGAGACTGT
TTAATAACCTGTCTTCCAGCCAATTGGTGGTGTGGAATCCCCTAGGAGCCTTCAGTCTGGGAGAAAC
AGAGCCAGACATAGACAGTTCAGCATCACAGAACAGAGAAGAGACCTGCAACTGTGAGAGTCCAGA
CAGGAAGCAGAGAAGGCGTCTTTCGCGAAAGGGCATTTCAGTGAGGCTTTGGAGTACGAATAGGAGCT
CAGCAGGCAGACGAATGAGGAATAAAGGTCAGAGAAGGTCAGAGCTGAGTGACGTTTGAATCCACCCC
GTTTATTGTAGAACTGGGGTTCAGAGGGCAGGTGCCTCAGAGTTGAGGCCACACAGTGAGGTCTGGTG
GGTGAAAGGACCCAGGAACGAGGCGTTCAGGAAAGCAGGTTGTCAGAGCTATGTGGAGTCTGTGGGTGG
CAGGGGCAGCCGCTCCAGCCTTTGAAGACTTTGAAAGCCAGAGATTCTTGGCGCAGGCTTGGAATTCCT
GGGAGCTCCTCCAAGTACCCAGGGGCATCAGAGCTGCCTGGGTGTTACATGGCCAGGGAACCCAGGTT
CAGGGTAGGACAGGCAAGACCAGATACCCAATGTGCAAAGTGAAAACACTGGGCTCCTGTAAACGAT
GAAGAATTCAAGACAGTGACAGCATTACGTCACCCCTGGGGACAGAGGTGAGCCTAAGGTGACACACGG
GGACTACTGTGCTTCCGGAGGCTCCCTGTGTCCTGGAGGAGAAAAGCATTAGAGGGGGCAGCTGGACAA
GCTCCCAACTGCAGAGTCCAGCCCTGGCTGGGGCAGGGCCCCGGCTGGGACTCAGCATTTCTGATAT
GCCTTAAGAATTCATTCTGTTTTGTACAATTATTTTTTAAAGTAAACGTGTGGAGAAAGAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001025604.3</u>
Locus ID:	27106
MW:	43.1