

Product datasheet for SC212791

LRRC10B (NM_001145077) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LRRC10B (NM_001145077) Human 3' UTR Clone
Symbol:	LRRC10B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001145077
Insert Size:	1171 bp
Insert Sequence:	>SC212791 3'UTR clone of NM_001145077 The sequence shown below is from the reference sequence of NM_001145077. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CTGGAAGCCGCTCCAGGACTGGGCACCTTGAGCCTATGTTCTGGCTGATGGGCTTCAGCCACTCTGAGAG
GAGGGGCTCTAGGAAGCTTTGGCTGGTTAGGCCTGCGGGACTGGTGGGTCCTACTTAGGTCAATGGAA
CGGCTACTTTTGGCGAGTTGCCGGGCACAGGCAGGGCTGGGGCTCATCTTCAGACACTCCTAGGTAGTA
AGAGGACTGACTATACCTGGCTGATCTGTGGGCAATTGTCAGACCAAGAAGTCTGGGCTGGGTCCAGT
CCCCGATGGAGCCACTGATAACAGAGGCCACAGCTTCCTCCTGGAACCAGGGTCTCTGAGGGCCACAGA
ATTGTACGTCTCCCTGAGTCTGGATGCAGGGCACCCAACACCCTCCCTCTTTCTGGGGTTTGGGTCTC
CTAGGCCACCCCTTTCTAGGACAGGAGCCGCTACCTCTGGGATGGACTGCGTGCCTGTGTGGGGCAGGA
CTGCCCTTGCGAGGAGGGCCTGGCCAGGAGCAGGGGCTTGGAGTGGGGCTGGTGGCCGGACAGAAGT
GGGAAATAGATGGGGGCCCAACACCGAAGCTGCCCTCACCCCTCCATCAGGGCTGATCGTACATACCT
CCCTCCTCGGCCCTGAGATTGGTATGGGGCACAGAAGGCAGGGGGTCACCCCAACCCTATTCTGTTC
TCTTCTCCAGCCTCCAGCATGTGCTGGAACCTTGCTGGGAAGAGAGAGCTCCCCCTCCACTCCCCA
CCCCATAGCCTGAGCCTTCCAGCTGTCAAGTTTGAAGAAGCCACCGTGCCTCCAGGGCCTCCTGCAC
CGCAGAGAGGGGCTCTGCTCTGGGAGGCTGCTCAGCGCTCTCCACGGCCTGCACTCTCCTTTCCC
TTAAATGGGCACAGCCAGGGTGGCTGGCATGGGACTGGCACTCTGTACCTCCACTCCGTCACCAAGGAC
AGCCCCGAGGGGGTGAGGAGAGGGGCTTGGCTCTGCAATCCCAGAGCCAGATCTTCTCATCCAAGCA
ATAACGGGACGAGGGAGCCAGGGCCAAGGACACTCAGGGACACTGTTGGGGCAGCGGCCACATGGTA
ATGTTTTCTTTAGAGAGAAGAAAAAACAACCTGGCTGTAACGTAATATATTTCTTGAGAAAGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001145077.2</u>
Locus ID:	390205
MW:	41.4