

Product datasheet for SC212812

MXRA5 (NM_015419) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MXRA5 (NM_015419) Human 3' UTR Clone
Symbol:	MXRA5
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015419
Insert Size:	1181 bp
Insert Sequence:	>SC212812 3'UTR clone of NM_015419 The sequence shown below is from the reference sequence of NM_015419. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TCCAAAACAACCTACATCCACGTCTTCTGAATGTGGATTCCAGAATGATTGCTTAGGAACTGACAACA
AAGCGGGGTTTGTAAAGGAAGCCAGGTTGGGAATAGGAGCTCTTAAATAATGTGCACAGTGCATGGT
GGCCTCTGGTGGGTTTCAAGTTGAGGTTGATCTTGATCTACAATTGTTGGGAAAAGGAAGCAATGCAGA
CAGGAGAAGGAGGCTCAGCCTTGCTGAGACACTTTCTTTGTGTTTACATCATGCCAGGGGCTTCATT
CAGGGTGTCTGTGCTCTGACTGCAATTTTCTTTTTTTTGCAATGCCACTCGACTGCCTTCATAAGCGT
CCATAGGATATCTGAGGAACATTCATCAAAAATAAGCCATAGACATGAACAACACCTCACTACCCATT
GAAGACGCATCACCTAGTTAACCTGCTGCAGTTTTTACATGATAGACTTTGTTCCAGATTGACAAGTCA
TCTTTCAGTTATTTCTCTGTCACTTCAAACTCCAGCTTGCCCAATAAGGATTTAGAACCAGAGTGAC
TGATATATATATATATTTAATTCAGAGTTACATACATACAGCTACCATTTTATATGAAAAAAGAAAAA
CATTTCTTCTGGAACCTCACTTTTTATATAATGTTTTATATATATTTTTCTTTCAAATCAGACGA
TGAGACTAGAAGGAGAAATACTTTCTGTCTTATTAATAAATAAATTATTGGTCTTTACAAGACTTGG
ATACATTACAGCAGACATGGAAATATAATTTTAAAAAATTTCTCTCCAACCTCCTCAAATTCAGTCAC
CACTGTTATATTACCTTCTCCAGGAACCTCCAGTGAGGGAAGGCTGCGATATTAGATTTCTTGTATGC
AAAGTTTTTGTGAAAGCTGTGCTCAGAGGAGGTGAGAGGAGAGGAAGGAGAAAACCTGCATCATACTT
TACAGAATTGAATCTAGAGTCTTCCCCGAAAAGCCAGAACTTCTCTGCAGTATCTGGCTTGTCATC
TGGTCTAAGGTGGCTGCTTCTTCCAGCCATGAGTCAGTTTGTGCCATGAATAATACAGACCTGTT
ATTTCCATGACTGCTTTACTGTATTTTAAAGTCAATATACTGTACATTTGATAATAAATAATATTCT
CCCAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015419.4</u>
Summary:	This gene encodes one of the matrix-remodelling associated proteins. This protein contains 7 leucine-rich repeats and 12 immunoglobulin-like C2-type domains related to perlecan. This gene has a pseudogene on chromosome Y. [provided by RefSeq, Mar 2010]
Locus ID:	25878
MW:	45.5