

Product datasheet for **SC212835**

PP14397 (JMJD8) (NM_001005920) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PP14397 (JMJD8) (NM_001005920) Human 3' UTR Clone
Symbol:	PP14397
Synonyms:	C16orf20; PP14397
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001005920
Insert Size:	1154 bp
Insert Sequence:	>SC212835 3'UTR clone of NM_001005920 The sequence shown below is from the reference sequence of NM_001005920. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGCGTCTTCATCTCCACCTTCTCGGCAGCAAAACAGCTGGCAGGACTGCCGGTCACACACCAGCAC
GTCCCACCTCGTGCTCAGGATTTTATTACACAGATAGTGGCGGCAATGGCCTCAGCCCAGCCCACCCT
CACCTGCTTTTCCAGCCACAAAGGGGGACGATCACGGCCAGCAAAAGCGATGCTGAGAGGGGAAACA
GTCCAGAGTCCAACAGCAGAACTTGGGGGAAGCGGTGGGGTGGCCAGGAACATAAACTATGTATAGGG
GCCGGGGGCTTCTGCCAGGGCTCCCCTGGACCAGGACGCCAGGTAGGGCAGGGAACCTCAGTAGTCCT
CCACCCAGCCATTCTCAGAGATGAATGCGTCAATAACCTCCTTCATAGCCAAGTTGGGGATGAGCTGTT
CCTGGGTGAGGGGGCTCCGGGTACGGGGTCAAAATGACCCACACGCTGCAGTGACAAGAAGGGCAGAG
GGCAGTCATGGGGCCAGGACCATGCCACTGGCCCTGCTCCCCAGCCGAGGCCCTACCTGCAGGTGC
TCCTCGATGTCCTTGGGTGCTAGGTGATGCCACTGGGCGTGATGCACGGCTCCCGCATCAGCTCAAAG
CTGATCTTGCCACACAGGTAGTCGGGGATGTCTCGCTTCTGTGGCACAGGGGCACACGGTCAGAGGCTG
AAAAGGGGCACTGCACGAGCACCTGCCAGCCATCGGCAGCAAGCGACACACACTACCTTCTCTCTC
ATCCACCTGAGAAAAAGCTCGTCCATGTCCGCCATGTACTTGTCTGTGAAGAGTTGAGTGCTGTGCT
TGGGGGAGACACCCACCTCCCTCCTCCATGGGGCAGACCCCAACACAAGCGGGGATGCTCCACGC
CACGTGCACACACAGACCCACATGTGGGTGGGGGGCACCCCTCACGTGCTTGGCCTCAATGCAGGCCT
GCTGGGGCCGACGTGGCTGTCGTCTCATCACCTCGTGGTTTCGCTGGCACTCTTCCAGCTCCCTGG
GGGTTGACCAGGAGCCGTCAGAGATGGACCTGGCCAGATGTCTGACCACACCCCAATCTCAGAGCTAA
CATCCACACTTCCCCACATTTCTGCTTGCCAGTAAAGCCTTCGATAAAC
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001005920.4</u>
Summary:	Functions as a positive regulator of TNF-induced NF-kappa-B signaling (PubMed:27671354). Regulates angiogenesis and cellular metabolism through interaction with PKM (PubMed:27199445).[UniProtKB/Swiss-Prot Function]
Locus ID:	339123
MW:	40.9