

Product datasheet for SC212836

DDTL (NM_001084393) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DDTL (NM_001084393) Human 3' UTR Clone
Symbol:	DDTL
Synonyms:	DDT; KB-226F1.2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001084393
Insert Size:	1180 bp
Insert Sequence:	>SC212836 3'UTR clone of NM_001084393 The sequence shown below is from the reference sequence of NM_001084393. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAGGAAGCTCTTTCATTTATTTTCATATGAGGATGAAGAAGAGGATTATGTGATCACAGGAATGTTGCA
TGCGGGATAATCCAAAGCTGGTTATCTCCAGGCCCTCACTCTGCCAAGAGATCTCTCTGGAAGAAGCAG
CCAGTTCACAGATGCCCTGGATCCCTCCGTGCCCAATCATAAAAAAGTCATGACCGTCCCTATCTTGCC
AATCTGCCAGGACTCCAAGGGGAAAAAGCGGATAAGTATCCTTCAGGAGACAGAGAAAAAGATATCATC
AGCTCCTTGGCTAATACCACATCTTGCAAGACCCCTGCCAGGTACTCCCACTGTGGGTACTCAGGACAG
CCTGCCTCAGTCCACCAGGCATTTTGCAAACCTGCTCATCCCAATAATGATTTTCCCAACCCCCAGCA
GGGCAGTGGGACACTCAGGTGTGGGAGGTAGCCGCACATCATGACCCAGCTAGGACAGACACACAATAC
AGTAGCCTCGGTGTGTGTGCCGTACACTCTATCATCTCCATCAGTTTGTGTACTGAGGGGTACCCTGCC
TCCTAATTGCAGGCTCCCTAAGGGCACAGTACCTCAGCCACCTGTTTTCTCAAAGGCTGAGCACAGGGT
GACCCCCACCCACCAACCCGCCACTAATGTCTGGGCCATAAGTGAACATGCCCTCTCAGAGGTAA
CAGCAAGTTTCCAGTGAGGAAAACCAAGTGGGACAGCCTGGTTGGGCGGGCGGGAGTATGAACAG
CCTGCTTCTCATCATAAAATTGGCATAAATGGAGATTATCTACCTCATGGGGATAGCATGCTTGTCTG
AGCAAAGGGAGAGTGGATAGGGATTGTAAGGATCAAATGGGATCATGAGTGTAAGTACCTGTAGAGGA
CCTAGCACATGGTAAGTGATATAACACTTGTGTAGAGTAACAGTTTTCCCTGAGTATCCGTGGGAGAT
TTGTTCTAGAATCCCCAGAGGAAACCAGAATCTACATTTGCTCAAGTCCCTTACATAAAATGGTGAAT
ATTTGCATGTAACCTACACATCCTCTCATATACTTTAAATAATCTCTGGGTACTTATAATACCTAA
TATAATGTACGTACTATGTAATAGTTGTATTGTTAGGGAATAATGACAAGGAATAAAGTCTATACGT
ATTTTCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001084393.2</u>
Summary:	May have lyase activity.[UniProtKB/Swiss-Prot Function]
Locus ID:	100037417
MW:	43.9