

Product datasheet for SC212853

KIF16B (NM_024704) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KIF16B (NM_024704) Human 3' UTR Clone
Symbol:	KIF16B
Synonyms:	C20orf23; KISC20ORF; SNX23
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024704
Insert Size:	1179 bp
Insert Sequence:	>SC212853 3'UTR clone of NM_024704 The sequence shown below is from the reference sequence of NM_024704. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TTTGACTACAGCAGCCACGGGACGGGGTAGAGCCAGGGGTGATGGAGGAACCACCACAGCAGTGCCTTC
TCGTCGAAGCGGGCTCCGATGCAGGGCAGCTCCCCATGCGAGGATCCGGGTCTGCCTCCTCCTGCTGA
AGACAGACATGCAGCAGCGGGCCGGGCCACCTCACGTTTCCATACCTAGTGCCTGAGTTTGGGGATGG
GATGCTCTGCCTGCTGATGTGGCCCTGACAGGCAGCCGTTACCGTTCCATTGCGGTTGAACGTGGCCTT
TTCCACAGTGCTTCTTCTCACTGCGCAGCAAAGTTCGTCCTGTGGCAAGATAGATGTGGTTGGGC
CATCGTGGGTTCCCTGAGCCCAGCCAGCCTGGGACCTCCCAAAGTGGGTGGCTTACCAGACCACCTTA
AATGACTTTTCATCTGGTTTCCTCTTTTACCAAAATATACTCGTATTTTTTATATTTCTTCCATGTGGCT
GGCTATATTCCAAGAAAAGCATTTTAAATTATTTTATTGATTTTTTTCTTTTTTCCCTCATTTGAAT
CAGAACTTTTATATAAAACCCAAACACTGATGTTTACACAGAATTTATATTCTGCAAAAGGATTTTT
TGATCCAATCATGACTGTAGTCTTCCATGCTTGACAAATTGGATGTAGACAACATTACTTAAACTTCT
ATAAATCCCTACAATTAGGATATTTATTTAACCTTGAATATTCAAGAACATTCTCCCAATCTAAATGG
CTACTGTGCATTCTTGAGCTTTTTCTGCTAAGCACAATGAACGCAAGCTAAATGCATATTTTTAAG
TATTATTACATTTTTTTGTACAGAATCTATTGGATCTTTGGCTGGAAAACCTAGAATTTATAGCAGTTT
ATTAATGATACCTTAAATTACTCAGGACTTAATGTAGCATTGCACTTCTGTGTACAGTAAAACCTGCTTT
GTTTTACTAAAGAGAAAAATGTGAGTGGAAAAATATGTATGTGGTATATACTCAAATGTATATAATTC
TATCTATAGATTATATATGTATACATTCTGTACAGTAGTTCCATCAAAATATGAATAATTCACACCA
ATTTTATTAATGTATTTGCTTTTTCAAATTTAAATTGAGCTGCTATCAATATTAATGAAGTTATGG
CATCTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024704.5</u>
Summary:	The protein encoded by this gene is a kinesin-like protein that may be involved in intracellular trafficking. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2010]
Locus ID:	55614
MW:	45.5