

Product datasheet for SC212879

RIBC1 (NM_144968) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RIBC1 (NM_144968) Human 3' UTR Clone
Symbol:	RIBC1
Synonyms:	2610028I09Rik
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144968
Insert Size:	1196 bp
Insert Sequence:	>SC212879 3'UTR clone of NM_144968 The sequence shown below is from the reference sequence of NM_144968. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGC
CAGCCCGGCCCTCCAGACTCAGAGACCTGAGGCCTAGTGGGGATCTGTCCTAAGCTGAGTACAGCAGCC
AGGGGACCAGAACTGTGTGGGGAGGTCCACTTATCCATAGGGTTGGTGTGAAATGTTAGGTTGTGCAG
ACCCCAACCTCCCTGACTTCCTTCAGTAACTGACAGTGTAGCATTGTTCTTGATTAGTTTGGCACTCCA
TTACGTGGAGGCAGCACAGTGCAGACAGACCTAATGGCTCTGCCACTTACTAGCTGGTTTGGCAAATG
AATTTTCCCTTCTGTGCTTCAGTTTCTTCATCTTTAAATGAGAATAATAATTTATACCACATAGAGTA
GCTACTATTATGATATTGTCTGTGGCCACATCCTTCTCTCCATTAGAAGAAGCCCTTCAAGGGCAAGG
ACCATGTCTTATCTACTTGCACGTGTTACCCCTTACATCAACCCTGTTACAGCACTCTGTTCTGTGGA
CCACTTATTGTAGTTGTAGCTAATTGTGAAAGTATTTGTTAATATCTGTCTTCTCCCTAAATGTGAG
CCTACTGAAGATGAACCATCTCTTTTATGTCTACCCATCTCCAGTTGCATAGTACAGTGTCTGACAGG
TCTTAGGAACAGCAAGGAGCTCAGTGTGGCTGGAGCAGAAAAAATGTGAGTACCTTTGCTCATGCCA
TGCCCTCAGCCTGGAATGCCCTTTGCTCTTTTAGCCTGTATACAAATGTATTCAAGGATCACCTCTT
CCTGGAAGCCTTCTCTTTTCCACCATCCAAACCCCTTAATACCCTGTGTACATCCCTGCCATTGCT
CTTGCCAGATTGTGGTCTATGTTTGATCTGTGTCCCACTAGACAGGGAGTTCTCTGAGGACAGGATT
TAATTCGATGGTATCATTAGAACCAGCACATGGCCCTTAGACTGAGATTGAGAGATATGAAATGATTT
ACCCAAAATCACACAGCTCAGAAGCAGCAGAACTGACATTGGGACCCAGGTCCATTGGAGGCATCTTTT
AGTTATGCGCTTGAGAAATCGTTGTTGTGTCACCTAGTTCCCTCCACATCCCTACAACCTCGGACACTAA
CCCCCTCTAAACATTCTGACACACACAAACACAAACTTGCTGAACAAATGGCTGGATGAGTGAACA
GGAATGGAGCAAGCCAGTACAGT
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_144968.4</u>
Locus ID:	158787
MW:	43.8