

Product datasheet for SC212900

SPATS2 (NM_023071) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SPATS2 (NM_023071) Human 3' UTR Clone
Symbol:	SPATS2
Synonyms:	Nbla00526; P59SCR; SCR59; SPATA10
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_023071
Insert Size:	1200 bp
Insert Sequence:	>SC212900 3'UTR clone of NM_023071

The sequence shown below is from the reference sequence of NM_023071. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ACCTCTCAGACTGAAGCCGTGAACCTCTGAGAGAAAATCCAGTTGGCCTCTCTCCTCTATCCACACAAT
TCAACTTGATACTGGACTTTAGGAACTTACAGTTAGATGTAATAACAAAAAGAAGTTTATGCGTATC
ACTTTTTGTGCCATTCTAAGTATTTTTGGTTTCTTGCTCCTTATTTCTCTTTACCATTTTTGGAGGG
GAAGCTATTTTTTCTTGATCTTTCCAGTGATATGGATTGAATCTGGTTGGTCATTTCCACCAGGA
ATCAGAGGCAGCTGTTACCAGGTTTCGGCCTTCCAGTTGATCCCTCCACTGTCCAGGCTGTCTCAGGAG
GAGGTGAATCAGAGCTAGTCTGTCACTTTCCAATCTAAGCCGAAGCCACAGGTGCTGACAGGTTCCC
TTCTAACAACTATTTGACCAAGAGAGGCAATCAGGGGTTGTATACTGCACTGGATTGCGCAGAGAAGG
GAAAAATTAATTAGTGAAAAGGAAAGAACTAGGAAACATTTGAGAACCTGCCTCTATCCAGAATGT
GCTGGAGATTTGACACTCAATCAGTGTTTAGTCTTCTGCTTGGCACCATAGCTTAACCTGCAGTTTCT
TCAAAATGCCAATGCCTTGTTTCTATTACCTTAGATTGCAACCAAGTCTAGGGAAGTCTATGAGAAA
GTAGCATTTAATTAAAGTTTAAAAAAGGTTGGCGTTGTGGCTCATGCCTGTAATCCCAGCACT
TTGGGAGGCTGAGGCGGTGGATCACTAGGTCAAGGAGTTCAAGACCAGCCTGGCCAACATGGTGAAACC
CTGTCTGTAATAAAATACAAAATTAGCTGAGCATGGTGGCGTGTGCCTGTAATCTCAGCTACTCAGG
AGGCTGAGGCAGGAGAATCGCTTGAACCCAGGAGGCGGAGGCTGCAGTGAGCTGAGATTGTGCCACTGC
ACTCCAGCCTGGGAGACAGAGCAAGACTCAGTCTCAAAAAAAAAAAAAAGCATTTTTCTGTTTTAAT
GTTTGGAGACTTTTTTTTTTGTCTCTCCCTTTCTGAACCTATTCTGATCCTGAACCAATTTATAG
CAATATAATTAGTGTTAATCTAAGGCATTACTCATCAAAAAATTGCTGCAGTCTTTACAGTTGAATAA
ATAAAAACAATGCATAAATATGCCCA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_023071.4</u>
Locus ID:	65244
MW:	44.7