

Product datasheet for **SC212901**

DEDD (NM_001039712) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: DEDD (NM_001039712) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: DEDD
Synonyms: CASP8IP1; DEDD1; DEFT; FLDED1; KE05
ACCN: NM_001039712
Insert Size: 1203 bp
Insert Sequence: >SC212901 3'UTR clone of NM_001039712
 The sequence shown below is from the reference sequence of NM_001039712. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGGAACCTGATGCTGCAAGCATTGCCCTGACCTATTCCCTCTTCTCACTTTGGGGACTGTTCCCATCAC
CCACCTCTGGAGCTTACACTGTTCTGGGGTTTGTCTCTACCCTTCCAACCAATCACACCCCTGCCTTT
TTTTTTTTTTTTTAAAGAAAAAGACAAAAGAAAGTGGAAAGTGGTATTCCCACCCCTCCCTGCACCCA
TGTGCCTGGGCTTCCCCTTTATTCCCTTTTCCATTTACCCCGTAATGTGTCTCTACAGCTACCTTACC
ACTGAGCCGTAAGACAAATGTATAGGAAGAAGCAAAGTCTACAGCACATAGTCTTTGTAAGGGATTGAT
GTGAACACTTTTTTTGGATGCACTAAGGAGTTATCAATACTTCTGGCTTTATGAGAGCTCTAAATTT
TGTCTAAAAAACCAAGGGCTGTGAGTAAGGGAGCTATGTGAAAGTGGGACTCTGAAGTGTATTTTGA
AAATTAATCGCCACCCTCTTCAAATTATAGAATTTTTTAAAAACAAGCTGTGGCCCTTCCACTCTCT
CCTGGCCTCTGGTGTCTCTCTCTGCGCCTTCTCCATTCCATGGCTTGAAACCTCTGCCTGATG
TGGCTCCTTCTTTTTCCATTTGTCAAACCTCTTCAAAGGAGGAACAGATAAGACTGGGTGAGCCTA
GTCATGCCTCCCAACTGTGGTGGTATATGTGTACACACATACACAGGGTGGATGGAGAAGAGGTTGCTG
ATTAATAACTACTCCCTTATAAAGGGGAGGGGGAGTGTGACACTTTCTTTCCATGTTCAAGTAAAA
TAAATAATGTACCCTGCAGCCCTTTCCCTTTGCTTTCTTCTGGCTTGGGCAAAGGGCATCATAGGTG
TAAGTGAAGTAATTCCTTTTTCTTCCCTCTCCACTCCTACGCCACTCCCATGCTTGGGAGAATG
GGGCTGGGACATGCACTGAGTGTGCACTTTTATTTAGGTAGGGAGGCATGTTGAATGAGCCAGGAGGT
CTAGAAGTGGAGCTTAGTCCAGCTGGTACAATACCAACTCCCCTTCTAGTTCCCAAAGGCATGTCCAG
ACACAGACTTTATGATTATTATTTTTCAATGCCAGTGTCTCAGCCCTCAGCAGAATTCAGTTTC
CATGAAATAAACAATGACTATATAAAATTC
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001039712.2
Summary:	This gene encodes a protein that contains a death effector domain (DED). DED is a protein-protein interaction domain shared by adaptors, regulators and executors of the programmed cell death pathway. Overexpression of this gene was shown to induce weak apoptosis. Upon stimulation, this protein was found to translocate from cytoplasm to nucleus and colocalize with UBTF, a basal factor required for RNA polymerase I transcription, in the nucleolus. At least three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	9191
MW:	45.3