

Product datasheet for **SC212902**

DEDD (NM_032998) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: DEDD (NM_032998) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: DEDD
Synonyms: CASP8IP1; DEDD1; DEFT; FLDED1; KE05
ACCN: NM_032998
Insert Size: 1203 bp
Insert Sequence: >SC212902 3'UTR clone of NM_032998

The sequence shown below is from the reference sequence of NM_032998. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGGAACCTGATGCTGCAAGCATTGCCCTGACCTATTCCCTCTTCTACTTTGGGACTGTTCCCATCAC
CCACCTCTGGAGCTTACACTGTTCTGGGTTTGTCTCTACCCTTCCAACCAATCACACCCCTGCCTTT
TTTTTTTTTTTTTAAAGAAAAAGACAAAAGAAAGTGGAAAGTGGTATTCCCACCCCTCCCTGCACCCA
TGTGCCTGGGCTTCCCCTTTATTTCCCTTTTCCATTTACCCCGTAATGTGTCTCTACAGCTACCTTACC
ACTGAGCCGTAAGACAAATGTATAGGAAGAAGCAAAGTCTACAGCACATAGTCTTTGTAAGGGATTGAT
GTGAACACTTTTTTTGGATGCACTAAGGAGTTATCAATACTTCTGGCTTTATGAGAGCTCTAAATTT
TGTCTAAAAAACCAAGGGCTGTGAGTAAGGGAGCTATGTGAAAGTGGGACTCTGAAGTGTATTTTGA
AAATTAATCGCCACCCTTTCCAAATTATAGAATTTTTTAAAAACAAGCTGTGGCCCTTTCCACTCTCT
CCTGGCCTCTGGTGTCTCTCTCTGCGCCTTTCTCCATTCCATGGCTTGAAACCTCTGCCTGATG
TGGCTCCTTCCTTTTTCCATTTGTCAAACCTCTTCAAAGGAGGAACAGATAAGACTGGGTGAGCCTA
GTCATGCCTCCCAACTGTGGTGGTATATGTGTACACACATACACAGGGTGGATGGAGAAGAGGTTGCTG
ATTAATAACTACTCCCCCTATAAAGGGGAGGGGGAGTGTGACACTTTCTTTCCATGTTCAAGTAAAA
TAAATAATGTACCCTGCAGCCCTTTTCCCTTTGCTTTCTTCTGGCTTGGGCAAAGGGCATCATAGGTG
TAAGTGAAGTAATTCCTTTTTCTTCCCTCTCCACTCCTACGCCACTCCCATGCTTGGGAGAATG
GGGCTGGGACATGCACTGAGTGTGCACTTTTATTTAGGTAGGGAGGCATGTTGAATGAGCCAGGAGGT
CTAGAAGTGGAGCTTAGTCCAGCTGGTACAATACCAACTCCCCTTCTAGTTCCCAAAGGCAGTGTCCAG
ACACAGACTTTATGATTATTATTTTTCAATGCCAGTGTCTCAGCCCTCAGCAGAATTCAGTTTC
CATGAAATAAACAATGACTATATAAAATTC
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_032998.3
Summary:	This gene encodes a protein that contains a death effector domain (DED). DED is a protein-protein interaction domain shared by adaptors, regulators and executors of the programmed cell death pathway. Overexpression of this gene was shown to induce weak apoptosis. Upon stimulation, this protein was found to translocate from cytoplasm to nucleus and colocalize with UBTF, a basal factor required for RNA polymerase I transcription, in the nucleolus. At least three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	9191
MW:	45.3