

Product datasheet for SC212925

LUC7L2 (NM_016019) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LUC7L2 (NM_016019) Human 3' UTR Clone
Symbol:	LUC7L2
Synonyms:	CGI-59; CGI-74; LUC7B2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_016019
Insert Size:	1144 bp
Insert Sequence:	>SC212925 3'UTR clone of NM_016019 The sequence shown below is from the reference sequence of NM_016019. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TCTGAAGAGCGCGAAGCAGGGGAGATCTAAGTAGCTGTGTACATTTCTTCAGTCCTTAAGCTTCCTACG
GAGTTACGTAATTTGTTAGTTACAGCTGTTACAGGTGACAGTGAGCAGATCCAGACACCAGATCCA
GCTAGGCTAGATGTACAGTATCTAACTTGATCTGAACCTGAACCTGTTTCCTTGATGATGCCTAAACT
ACATCCATAGTTTCTGGTGAACCTGTAATACAGTTCTGAAAGTACAGTTTATATAATAAGATGCTGAT
CTCTTTATTCTTTCAAGTAAGAGTGCTAGTGAACAAATTGTGTTACTTGCCTTGGGATTTTTTGAACGT
TTGTAAATGCTGTCTTCTAGTCCAAACAGCTGCAGCTTTGGGCATTTTTCTTTTAATTATTCTTCC
TCTGACTTTGTATCCCTTAATACCTACACTCTCCAATTGTAAGAGAAAGGGGCGAGGAAGCAATATAG
CTTCCATTCTAAGGCTGTATTCCCGTTATGAATTACTAGCTGATTACAGTTACAGCATTGATCCTGGA
ATGTGTGCTGGAGAAATTTAAATACTGGGGTTTTTTGTTAATGGTGCCTATTTAGAGTTGGAAGTTG
AACAGCTGTTGCATTACATACTTTTGTCTTTTATTGAAATTTTGAAATCAAACGCTTGATTTTTCTG
TTCTGTTGAATTGCTATGTTCAAGATGTTCTAGGGGTGGGGGCGAGGACTCTTTTCGTAATAAGCACT
TGTTTTATTTTGTGTGTGGAGTATAAAGGCTACACCCTTATTGTAAAAAATAATAATAATAAATG
AAAGAAACAATCACCACCACCATTATTAACCTGTAACCTGTCTAGTAAATTGTCACTAGAACTATTGTC
TGTGGGCAATTTCTCTATTCTGTACGTACATTAACAGTGCCTTACTGTGCAAGGCACAAAGGACATAA
ATACGTAATGCAAGATTCCAGATGAGCTGTTTATGTGGGCCCTTGTGACTATATTCCAGCCACTT
CAGGGTTGTTTGAATGACCGTTATAGAGAAGGGCTCGACCTGCAGAAGAACTGGGGTGTTTTTATTC
TCATTCAACAAACCCAAATAAAACAGTATATGTAACCAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_016019.5</u>
Summary:	This gene encodes a protein that contains a C2H2-type zinc finger, coiled-coil region and arginine, serine-rich (RS) domain. A similar protein in mouse interacts with sodium channel modifier 1, and the encoded protein may be involved in the recognition of non-consensus splice donor sites in association with the U1 snRNP spliceosomal subunit. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2011]
Locus ID:	51631
MW:	42.7