

Product datasheet for SC212927

C17orf108 (LYRM9) (NM_001076680) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C17orf108 (LYRM9) (NM_001076680) Human 3' UTR Clone
Symbol:	C17orf108
Synonyms:	C17orf108; HSD24
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001076680
Insert Size:	1145 bp
Insert Sequence:	>SC212927 3'UTR clone of NM_001076680 The sequence shown below is from the reference sequence of NM_001076680. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATCATGAACAAATATAAAAAACAAACAGAGCTCCGAGGGTGAGCCTTCTGGCCTCTGGAAGTTGA
GAGGCCCTCTTCTGGAGCTCAGCAGCAGCAGCTGGTTTGAAGAGCCATTGGCCTTGTGGCTGAAAG
CCAGAAATCAGGGGAGAAGCTGACCTTGCTCCAGAGGCTCCTTTGCCGCGGCCATGTTGCTAGCATCC
TTTATGTTTGCACTTCCAGCCACTTACGATGTGAAATTTGGGAGAATGTAGTATTTTTTGCCTCTTT
TGTGTGTTTTTATTTTTGTGTTTTGTTCTTTACTTTCTCTGAAAGTGATTCAATTTATTCTGGAGA
TCTCTTTATCATTCTCTGAGTCCAGAGCAGATGTCCCTCAGGTAGCAACCTGTACAGAAGAAAGCAC
GTGGGTTTTAGAACAGACAGATCTGGATTTGAATCCCAGATGCAATAAGATGAAGGTGAAGTGCCTCA
CACAGGATGTGGCAGATGACAGTGCCTGAAGGGCTCCCTTCAAGGAAGCCGTCCTGCCCCGCCCTCC
TCTGCTCCTGCAGCGCCTGCACAGCTCTGCTCCCGCATGCACCGCATCGGGTCTCAGTTTTTGTGTT
GCGAGCTTGCCCTGTCCCTGCTTGTAGCCCTGAGCCCTCACAGTTGTACACACCCAGCAGTCTGCCA
AGGAGCTGCTTTTCTGACTACTGTCTCTCCATGCCTAGGGAAGTCAAAGCCCTGGGCTTGCTTCAG
GGAGTTGACAGTGTCTTGAAGTGGTGAAGTGCCACCATGAACCAAGGAGAGACGCCAAAGCAGCT
GCCATGGACTCTGGCCTGTCCAAGCCTCTGTGCACCTGCAGTCTCTTTCCACTGGCACTTCCCTGCTTT
GCTACCATCGTGAAGAATGGAAGGGAGGGGAGAAGGTGGGTGGGCTCAGGGAGAACGGCCAGAGTGG
CCAGAAGATGCAGAGAGAAGCCCACTGGGATCTGAGTATTTCCCAAGGGCACTCAGCCACAGGATGA
ACTAGTGTGTTGTATCCTTGCAATTATGTGACTTAAAGTGACTCAGAGATATCTATGTATGCGAGT
GTACTCCGCTGGCAATGAAGTAAATATATCAATTTGTGGAA
ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001076680.3</u>
Locus ID:	201229
MW:	41.9