

Product datasheet for SC212944

SIPA1L2 (NM_020808) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SIPA1L2 (NM_020808) Human 3' UTR Clone
Symbol:	SIPA1L2
Synonyms:	SPAL2; SPAR2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020808
Insert Size:	1192 bp
Insert Sequence:	>SC212944 3'UTR clone of NM_020808 The sequence shown below is from the reference sequence of NM_020808. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TTTTTCACCACCATCGACAAAAATCTTAGCCAATCCGCACCTCATCAAGGGACCACTGGGAAATGCCC
CTTGTCCTTTGAAGTCACAAACATGTGGTTTTCTGTGTGCTTTCAACCAATCGTAGATGTTTTGCT
GTTCCATTCTGTGTAGCACCATTACCACAGCAGGATAGGGAGCCTCGACTCTTTCTCGGTAACCACGG
CAGAGAGCAGCGCCGATGTGAAGATGAATGAATGTAACCTCTGGTGTGAAGATGAATGTAACCCCTGGC
GACAGTTGAGACCTTTTTCTTTAGACTCTGCTAAACAGTGCTCTGGCTTGGGCTTACCTCAAGAGGG
AAGATAGTTGAGTTTTATTTCTGTATATCAGGTGACCTGGTAGAGATGTAAGCAATTTACCATAGTT
TGGGCTTTAGTATTGTAAATAAACATGAGAACAAATAATCAGACATATACTTTAATGTTAAAGGTGCT
CTATTTTTTTGGATGTACAGTAGTTTTATTTCCACAGCCACATTACCATAGCAATAAGAAAGGAGGCAT
AGTACATAGTTGGAAAAGCTTTGTGGGGGATTAAAAAAAAAAAAAGCACTGTTGTGTTAACTAG
TTCAGATGCAGTTACCTTAGAGACTTATTTATTTGCAGGAACAAATGGTGCCTGAATATTAACAGTGTT
CTGATTAACCAAAAAAGATACATATGCCTTGTAATGGCTACCGAGTGGTCAGTAGTCACCTCA
ACTCTTAGTTCACTTTTGTATAGTTGCTCTGCTGGAAAGAAATGAGAGTGAATCTGCTTACTCACTAGA
ACTTCCCTGTGTGCTGTGAGCCAGCGGAACCACTTGTAACAATGCCAGATTGTTTATCTTTGTACAGAA
GCTTTGATGAAGTGCTTGTATTTAACCCCTTATTTAAGCTTATTTAACCTTAAATTGTTAATTTTAT
AAAATTTGGTTTGGCCTGCACTACGATAGGGATGGAGGTAGCTGCAGGCTCAGAAGAGACTGAGCTTG
CACAGATCAGACCGAGAAGCAGGGTGAGAGATTCTAACGACTGGATGCTGCTAGTAACACATTGTTTGT
ATTGCTTTACATTTTAACTGTTAGATTTGAGATGAACATACATTTTGTCTTTTAAATAAATGTTTAA
AAGAAGTCCACATAAGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_020808.5</u>
Summary:	This gene encodes a member of the signal-induced proliferation-associated 1 like family. Members of this family contain a GTPase activating domain, a PDZ domain and a C-terminal coiled-coil domain with a leucine zipper. A similar protein in rat acts as a GTPases for the small GTPase Rap. [provided by RefSeq, Sep 2015]
Locus ID:	57568
MW:	45.8