

Product datasheet for **SC212950**

ABI3BP (NM_015429) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ABI3BP (NM_015429) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ABI3BP
Synonyms: NESHBP; TARSH
ACCN: NM_015429
Insert Size: 1190 bp
Insert Sequence: >SC212950 3'UTR clone of NM_015429
The sequence shown below is from the reference sequence of NM_015429. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGTGGGACTACAATTCTGGAAAATGGTAGATGCTGCACAAAGTTACCTTCTGTTTCATCATTGCAAAC
AAAAATCATTGAAAATACTATGCCGCATTCATTTAAAGCTATTTTGTACTATGTATAAAAGTCTACA
ATCTAATTAATAGCAATACTAGATGTTTTATTATTAGAAAAGATTGCTGAGAGTATTATCAGGTTTTAC
AAAGTCATTTTAAAGAAAGCAAGATACTGATGTTAACAGAATAACATTTTTGGGGAAGCTGGCTCCCTAT
TCATGGTATTTTAAAGAGATCATTGTATATTATTTATCACACTGTTGTAATGATGTTTTGAGATACTTT
TATAACAAAATTAACATCAAAAAGGTATATACTTTTTAAAAAAATTTACTTTTTATTGATGTGACTCT
TCCTATTGATGAGTTAATCCATAAATCTACTTAGTTTAACTTATTGGATCAAATTATCTTCAGCAT
GTATATCTGGGAAAAAAGTCCGAATTTTACATTTTATTTAACTTCAATTTTTTATTTAACT
TCAATTTTTTAGCAACAGCTGAATAGCTTTGCGGAGGAGTTTAAAGTTACACATTCATGCTAATATAC
ATTTCTTTAAACATCCACAAATTTCTTAAAAAGATTGAATCAGTAAATTTCAATTCAGCTAAAAATGGA
GTCTAATATATTGTTTCAAAGATACATTTTTACCACCATAAATGTTACAATATCTGAATATGCTTTG
TCAAACACTACCTTTATGCAATCGTCTTCATATTGTTTTATGATTCTAATCAAGCTGTATGTAGAGAC
TGAATGTGAAGTCAAGTCTGAGCACAAAAGATAATGCACGATGAGATTGCCACCTTTTATAGGATA
TTTACTATGTATTTATACGTTAAGACCTCTATGAATGAATGTATCAGAGAATGCTTTGTAAGTGTTTA
ATTCAATCTGTAATAAAATCTAACTAACTAATCATTATTTCTATTTAAAAAGGATTGTCCTTTAGG
CGGGGAATGGGAATCCTTGCTGCACTGTTGCAGTCATTCTGAAAGGACCTTTCCCTGTACTTACCTTTC
AACATGCTTCAATCTTATCAACGCTACATTTTGTATTTTTCAAACAAGTATAAATCTGCAATAAAGAG
ATGTAGTTTTTTTTTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015429.4</u>
Locus ID:	25890
MW:	47.2