

Product datasheet for SC212956

MYH15 (NM_014981) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: MYH15 (NM_014981) Human 3' UTR Clone
Symbol: MYH15
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_014981
Insert Size: 1207 bp
Insert Sequence: >SC212956 3'UTR clone of NM_014981
 The sequence shown below is from the reference sequence of NM_014981. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GAGTTTGGGAAAAAGGTTCAAGAAGAA TAGCATCCCCTGCTTTGAAAGGACAACAGCTGGAGAAGTACA
AGGAAGGTGCTGTTTCATGGCCAAAACTTAGGTTGCATGGAAACATTTTTAAAAACATGTTTAAATTG
CTTTTCACACCATATAACAAGGCAATTAGAAAAATAATTAAGGGAATATCATTGCTCCACAGTTAA
TGGGGATTTTTGATCCTCAAATGCAAGTAACTACCTTCTAATGCTTCACATGACAGATTAATAAAT
GGAAGAACCCTTTCAATTCTGATGTTTAAAAAATGAATAAATACTTGATCCTTTGTCCATATTTCTCT
TAATGGGTAGGACTCATAGAAGATGTCCTTAGACCAGTCACGCTTCATGGGGACTAGGGCATGTTGGTG
AATGGTTTTTACTAAAGTTAGGCAACTTTGGCTTGATTACCCCTAAATCTATGAATGTATATTGTGAG
CCAGCAGTGGTAGAATGGAAGTCATAAGATCTCCTCAGCATTGTGATATAAATATAGCTATAGTTAGGC
AATTTGAACATGTAGGCAAACTCTCCTAATCAACACACATGTAGGCTATATGCTGGTACATGCTTTAA
ACATGGAGGTAACCCACACGAGACATTCAGTGACAGGTAATATGCTGGATTTGTGCAGTACCGATTTG
GCTCAGCTGGAGAACATGTCCCGAATCCTCTCTCATGGTCTGGATTGGAGTAGGTCATAAAGA
AATTTGCATTAATTTGATCAGCAGCTATTTTATACTCATAAGGTCAGTGTACAGACCCAAGCATGGTGA
CAGCTTAAAAATATGACTCCAGGCCAAAAAGGGGAGCTAGAAGAGACCAGAGACAGCTCCCTGGACCCA
GAGCTCTTCCAGCTCCTGCCAGCCTCCTCTTACGCTTTGCAAAGTACTGGCCAGGTGTGTGTGCAGCT
CCATGGCAACCAGCATCAGCTTTTCTGAGATCACACACAGCATTGCAGTGGAGGCCGTGAGACAGACA
TGGGTTCTGTTTGTCTTATGGACTTCCCTTATCCTTGCTGCATTACAGTCGACCTACAGTGACTTC
AGGCCAGAACCAGATGCAGAGGGAACAGCCTGGCCTAGACTTCTCCACCAGCACCCACAATTGTGTA
GGCTGAAGCTCTATAATAAATCTTTATTCTGTCA
ACGCGT AAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014981.3</u>
Summary:	Muscle contraction.[UniProtKB/Swiss-Prot Function]
Locus ID:	22989
MW:	46.3