

Product datasheet for SC212965

UBAC2 (NM_001144072) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	UBAC2 (NM_001144072) Human 3' UTR Clone
Symbol:	UBAC2
Synonyms:	PHGDHL1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001144072
Insert Size:	1199 bp
Insert Sequence:	>SC212965 3'UTR clone of NM_001144072 The sequence shown below is from the reference sequence of NM_001144072. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GTCGCCACCAACTTCTGCTGCAGCACTGATGTCAGGCCAACACTGGGACGGACCGGCAGCCGAG
TGACAGTGCCTGGTCCCCACCATCAGATCAGCCCGGGACCGAGCATCTCTGGTCTGATGTTCTTG
GGAAGAGGGAGGTTCCACCGCACCCCTGCCCTCAACCGCAAGACTGTTGCCGTTTTAGTGTGGAGATAA
GTTTGCCATTACATTAGCATGTATTTCTATCTATATTTTTATTGGGCATTTTCCCTAGGTTGGAGAG
TCAGCACTCGTTTTGAATGTGTTTAAAAATGCATTAAATGGAAGATTTCTGCAGGCAGTTGAATGGCAC
TCCAGATGGGGAATTGCTGTAACCTCTTACTGTAACATGTCATCTCCTGCGTCGTGATGGGGAGAGGG
TAATGTTACTTCAAAAGGACATGTCAGATCCTTCTTCATGGACTTTTTAGTTACTGTTTTTCTCTC
AACTTGTTCGAATCTCCTGGGAGTGAGGGAGAAACAGGGAGCTGAATCCTCCCCAAGCTGTCCA
GGCCAGAGGACTCTGCAGTACCTTCTCTACATCTAGTAACAAAGAAATGGTGATAACCATGCACTGGT
CAAGGTTCTGGAGTTCTCCATGAACTTGGGTTAATTTTGCTCAGAGTATCCAGAGTTAGCCACTAGGC
TGCGGGTGAAATGGGATGGAGAAGAACAACAGCAGGCTTCTGGAGCCACATGGGCTGACTAGGGCACT
CTGTGGCTGGCCTGGCATGGGCTCAGCCAGGAAGAGGAGAAACGATCCCTTGCTGCCCTCCCTGTG
GCAGGGCTAACTGCCTGGCCCTCCTGGCTCGCAGCCAGCCAGCCCTGCGCAGCAGGTTCTCCTCAGGG
CTTGGGCTTCAACCTGTGGCGACAGGAGGCAGGGCAGACTGTGGAGGACAGGATGCAGGTGAGGGAGA
GGGAAGGCAGGGGTGGACCGCCATGAGCATGAAAAGACCCGAAGCAAGTTGACTCTTGCAATGTGAAC
TGTTATGTTCTGCAAAATGAGCAACGATGTATCAAATTGATGCAATTTAGATGTTGATACTTACAATA
AAGTTTTAATGTGTTTTACTCTTCAGTATTTTCTATCAGACTGTACAAATAAAGCCTGTTCCAGG
CTATACATATTTACAGTCAACCAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001144072.2
Summary:	Restricts trafficking of FAF2 from the endoplasmic reticulum to lipid droplets. [UniProtKB/Swiss-Prot Function]
Locus ID:	337867
MW:	44.2