

Product datasheet for **SC212968**

LOC100132288 (NM_001033515) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LOC100132288 (NM_001033515) Human 3' UTR Clone
Symbol:	LOC100132288
Synonyms:	FLJ00219; FLJ35473; FLJ39633; MGC90442
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001033515
Insert Size:	1208 bp
Insert Sequence:	>SC212968 3'UTR clone of NM_001033515 The sequence shown below is from the reference sequence of NM_001033515. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CCTCTGCACATAGCCAGACCCGGCTGTAGCTGCTCTCACACCGGCCCAACATGCAGCTGTGCCGCGAT
GCAGCCCAGTTCAAGTACTGCCTGGGCTCTGAGGCAGTCCCAGGTGCCCCATCCACCTCCTCGCTGC
GACCTCTCCAATAAACTCAATCTCCACACACATCCCCACAGGTGCTGAGAGGGGAGGGAGAC
TCTGCGGCCACACACAACATGCAAGCACGCTTAAACACACAGATATAGACATGACACTGGGCAGGTCA
TGGGTGCGGGGCTCCCCGTGCCAGGGACCTGTGCCACGTGGTCTACCCCACTTGTGTACCTCAT
CTTCCCGGAGGAAGCTCCTTCAGACAAGATGAGCTCCTGGAGGCCACGGGGCCAGAGAGGCCAGAGAG
GGCCATTCCCTGCTGCCCTGAGGGTGTGTCCCCACCTCTGGTGGAAAGTATGACCTACGCTGGGCTA
CAGGCATAGATGGGCGTCTTTTGTTCCTTGTAAACGGATGCCTGGTGGCCCTTCTGAGCCTAAG
TACCTTTCTTCTTCCCTGCCCTCAGTTCCTTCAAGGAAGTGGTGCCTCCTTCCGTCAAATATCAGG
AATGGAAGAAGTGAAGCTGGTGGCTGAGGCCAGTCTGTCTGTGAGGTGAGATTGCACGCGGGT
CACTCATTCCCACTCAGGCCGCGCTCCCTGGGAGGACCTCGTAATTTATGCTCTGCCCTAGAGCCGC
TTTCTGGTAGGATCCTTTCATGGCGTCTTCCAGGACCTCCTGCCAGTCCCATTGGGAGTTTTGGGT
CCCTCCACCCAGTAGTGCATACTTCTGCCCTCCCCAGCCCCAGCTCACACCCCAACACCCCAAA
GCTGGTGAAGTGAAGTGGAGGAGTGTACAAGTCCATCACAGCGCTGCGGGAGAAGCTTCTACAAGCGGA
GCAGTCCCTGCGCAACCTCAAGGACATCCACATGAGCCTGGAGAAGGACGTCACTGCCATGACCAACAG
TCTCTTCATCGACCGCCAGAAGTGCATGGCCATCGTACTTGTACCCCACTTCTGCAGCTGGCTGG
CTACAGTGAGCAGCGGCACGGTCTTCCCCCAATCCCCCAATAAACAGTGCCTTAGCTTTAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001033515.1</u>
Locus ID:	100132288
MW:	43.3