

Product datasheet for SC212971

ZDHHC5 (NM_015457) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZDHHC5 (NM_015457) Human 3' UTR Clone
Symbol:	ZDHHC5
Synonyms:	DHHC5; ZNF375
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015457
Insert Size:	1181 bp
Insert Sequence:	>SC212971 3'UTR clone of NM_015457 The sequence shown below is from the reference sequence of NM_015457. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGTGGTACCACCTATGAGATTTTCGGTGTGAGCCTTCGGCACCTCCCCTCCCCAACGCCTCTGCGCTAC
ACCAAAGGGCCCCAGGTGGCCACCTTCTTCCCTCAAGGGGCTCCCCTCCCGTGCATGGACATTTTTTA
AACCACCGATTCCAAGAGGATGAGGAGTGTCTTCTAAAATGCAGTAGGCTTGGGGAGTCGGAGAGTTGG
GGCCCTGAGACTGGGGTAGCAACCCCCCTTTATCTTTTAAGACCTTCCCTTCTTGATCCCTGGACC
AGACTCAGTGGACATTTGTGCAATTGCTCGCCCTGGAGGGAACCAGATCATTTTTAAACCAGAAATAAT
TTTTTTTATTATTGTTACGGATTCTATTTTTTCTCTCTGCGTTACCAGGTGTGTGTGTACATATAA
TATATATATATATATTATAAATATCAAGAAATTATATATCTATCCTGGGATGGGAAAATGAGGGAG
GGATACATATACGGAGGGGGATCTTACTCTTCCCATTCCTCAGACCAGCAGGAAAAGAGGGGAGACGTC
AGTCTTTTTCTGTGGTTCCTCTCATTGTCCAGTTACTAACTACGAAAATAGCATCTCTGCTGGT
GCTAAGTGTGATTAGGAAGAAGCCTGGGGAGAGGTGAGTCTGGAATTTTGGTCACAAGAGGGAAGGACT
TGGAGAGGAGAATTAGTTTTCTAGGCTCATTGGCATTAGTTTCCCTAGGAAAGGGTCAAACTTCAA
GACACTGGTGGTGGTGGGAGATCAGGAAAATAACTTGGCCTAGCTCAAACAATATTGGATAATCCCCTC
CTTGGGGGAGAGGATTAGAGTGTGCTCCTACTGGCCCCCTTGGAGCCTCCCCTAGCTTACACAGTTAAC
TTGATTTTAAATCCAAGGCCAGGAGAGAAGAATCCAAAAGCAATATTTTTTATCACATGCCAAAAC
GGGGGATAGAGAGAAGGAGTGGCAGGCCTAGGCCCTCCGATTGTCCCTTGGGGTTACCCCTCAGCCC
ACCTCACTATGGTGTGGTAGAGGGGATACCTGGGTTCTAACCTCTAAATAGGGGAGATCCCAGCCTC
CACAAAGAGGCCCTTTATTTTTTATTCTGATTAGCCATTTTAAACCAACGAGGAATAAAAAGAAATCC
TGATCTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015457.3</u>
Summary:	Palmitoyl acyltransferase for the G-protein coupled receptor SSTR5. Also palmitoylates FLOT2 (By similarity).[UniProtKB/Swiss-Prot Function]
Locus ID:	25921
MW:	43.9