

Product datasheet for SC212999

SNX18 (NM_052870) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SNX18 (NM_052870) Human 3' UTR Clone
Symbol:	SNX18
Synonyms:	SH3PX2; SH3PXD3B; SNAG1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_052870
Insert Size:	1237 bp
Insert Sequence:	>SC212999 3'UTR clone of NM_052870

The sequence shown below is from the reference sequence of NM_052870. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GTACCTTTGATGACAGTGCTATCATTTAGAGTTTGAATAGTTGAGTAATGAGCTGCTCAAATCTGTAA
TCTTGAGACGAGGGTAGAATTAGAAGGGGACCTAAACCATGTTAGGTATCTGTTTGAACACTGCACTGT
CACTTAGTGGTGAGTTTCAAAAAGTGAGGAAATAATGGAGAAGAACGGATGTCTGTAATAGAAATCTG
CTTTCACCCTTCCACACCTCTCAATCTCTGCTTTGGTAGAGGAATTAAGATCAATTCTCTGAGTTTAAG
AAATTCCTGCAAAGAATTCTTGCCTGCAACCTTTACGAAGTGAAGTAATTCAGGAAAAGCGTACTCAT
TTTGAGTTTGAATCTATTAATGTACTGAAATTGAGAACGTTGTTTAGGAAAATTGTTTAGGGAAACTA
CATAGAATGAATATAATCTGAGTATTGGATCTGAGAGCAGCGTTGCTGGTATCCCAGTGAGTTAGCTGT
CCAACATCAGAGTTCATGAGTTATTAATAGTTGTTAGTGTCTCTGACCTTTGCTAAAGCACCAGAATT
TGCAAATAAAAAATATATAGCAGTAGTGCTTGTCTTAACATTTCTTACTAATGTGCATTCTCATATA
ATGCCTTGATTGAGAGAAATGGGGAGGTGTGGGAGTTAGGGAGGCAGTCTTGCCTGTGAAAGTCCAGTGT
AGGGGCGGCTGGGCAGGAAACCTACTCCTGAGCAGGAAACCTACTGTCCCTGTTTCTCAGCGGGTC
TGGCCTTTACTGGTGAAGTGCTCACTGCCTTTCTGGTGGGCCAGTCCGTTTATGAGCAACCCTGATCTC
TGGCATTGCTTTTCTTGTATTGAACCTGAAGTTTGCCTCCCTTGATTGTCTTTCTGTTCCCTGGGACCGC
CTTCAGACATTTGGAGACAGCAGTCTCGTTCCTTCTCTGCTTCTCTCCTTCCCCTTCCAAATGCG
TGACTGGTCTTCTGACTCCTATGGGGGAAGGGGCTTGGTGGGAGAATGATCAAACGTTTTTCACTTGT
TGCTGCCGCTGAAAGGGAGTACCCTGGACATCAGTACTAATTTTAGTAGACAGTAGTAACTTTCTTGGC
AATTTGCGATTTGTAATGCAGTCTAATATTGAATTCGTTTTTGTGCAATCATATTGCTGTTGGCCTG
TATTAACCTGTATTTACATTACCTATGCCTTTTTTAAACATAAATTGTTATGAATCTTTGA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_052870.4</u>
Summary:	This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members, but contains a SH3 domain. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	112574
MW:	46.9