

## Product datasheet for **SC213028**

### TAF5L (NM\_014409) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** TAF5L (NM\_014409) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** TAF5L  
**Synonyms:** PAF65B  
**ACCN:** NM\_014409  
**Insert Size:** 1193 bp  
**Insert Sequence:** >SC213028 3'UTR clone of NM\_014409  
The sequence shown below is from the reference sequence of NM\_014409. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGAATTACACAAGAAAATCAGGAACATTAATTTTTTAAGTTTGTGGAACGGACTGGGGTAACAGTGGAA
AAGCCTCCAGTTGCAATCCTTCAATGAGCTGAGATTGAATCATTGACTGTGAAAGAGCCCCCTCCTCCC
CTCTCTCCTGGGCCCTCATTCCGGGCAGGTCCCCCTTCACTCAGCCCCCAAGTGTGAAAGCCAAGGTT
GGGGAGGGGAACAGGAGAGAGGCAATGATGAGTGATCAAGAATTAGAGTGCCTGGAGGTCACTCACA
GTCCTCATGGAGACTCTGCCGTGACCTCGCGCATGCGAGCACCTTCTGTCTCTCATTTCACCATCTG
ACTGGGGCTAGGTTTGGTTAGAGGGCGGCTGTAGAGCACCTCAGGAGAAGTGGCGCCGAGGGGTAC
GCATGGTCCAGCCCTCAGGTGGGAGGGAAGTCTTGTGGCAGAGGGATGTGGGGCGGGAATGAGGAAA
GGGGAGGGCATGATCCTGGCAAATAGTCTTTCCTTTCCGTAATTACAGAGAACCAGGATTTTTTTTAA
TTATTGTTATTATAATTATAATTATTATTATTATTGACATGAGGAAACCTAGCACTGGCAGAGGGCACC
TTCTCTGCTTTTAACTTTTCAGGTTTGACTCCTGGCACTGGTTTTGATACTTGGAAATTTGAGGTTTCAG
GGAATTCAGAGAATTTGGTAGCACAGTGTATTGGGAGAAAAGGGAATTTCCGGGTGACAGATTGATACC
CTGGCTGACATGTTTCTGGAAGTGGGTGGGTGGTTTTGATTGAGAACAATTTGGGCTCCTGCAGTGATG
GTGGGAAAGAAAATTTCTTTCCAAAGCTGCACTTCTTAATAATTCTCTCCCCTCAAAGAAAAGATTA
TGATAAAGAATAATTTGGAGTAATAAGTAACAAGATAGAATTAGAAAAATGTTGAATTTCTGAATGACC
ATTCAGAGGTGTGGCTTTTTCCACGTCTTTTTCTCCTGCAATGAAGTGTCCCACGTCTGCACTGCAA
TACCAGTCTAATGAACACTTTGTGAAACAAATCGTAAGTATGTACGCACGCATATCAAATATGCTCTG
AGGATAAAGCAAGTTCTAATTTGTGAAAAAGGATTAATATTTCTGTTTTCTGAAAAGAGTTATTTTG
TATTTTGTTTTCTATTAATA
ACGCGTAAAGCGCGCGGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_014409.4</a>
<b>Summary:</b>	The product of this gene belongs to the WD-repeat TAF5 family of proteins. This gene encodes a protein that is a component of the PCAF histone acetylase complex. The PCAF histone acetylase complex, which is composed of more than 20 polypeptides some of which are TAFs, is required for myogenic transcription and differentiation. TAFs may participate in basal transcription, serve as coactivators, function in promoter recognition or modify general transcription factors to facilitate complex assembly and transcription initiation. The encoded protein is structurally similar to one of the histone-like TAFs, TAF5. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]
<b>Locus ID:</b>	27097
<b>MW:</b>	44.5