

Product datasheet for **SC213034**

ATXN7L1 (NM_152749) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ATXN7L1 (NM_152749) Human 3' UTR Clone
Symbol:	ATXN7L1
Synonyms:	ATXN7L4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152749
Insert Size:	1218 bp
Insert Sequence:	>SC213034 3'UTR clone of NM_152749 The sequence shown below is from the reference sequence of NM_152749. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGCCAAGCCATAGAGAAGCATCAGGTGGAGAATGGAAAACGCAGAAGAGACGTACAACCTTCTGAAAGA
TCTCAGAGGACTGGCATTGTCTGGAGACGGCTTCTGGAAGAGGAGACTCCAGCAGAGCCCTGAAGTGT
GGACAAGCTTCCATTCAGTCAAATTCGTGAGCAGCTCTGGTGTGAGGATCTAGTAAGGAAATGGGCT
TGAGCAGAGCAGGGGTACCTGTGGGTTCTTCTTTCATTTTACCACCACAAAGCTGGTCTAAGCTG
CCATCAATTCACCTGGACAACCTGCAAAAGCCTTCCAGCTGTTCTGCTTCCACTCTTGCTCCTCTAC
AATTCATTCTGCACGCAGAAAGTCAGAATCAAAATCAAAGTTCATCATGTCACTTCTTGACTTGTGGCT
TCCCTTTGCTCTTAAGCTAAAAGTAAAAATCATTATCATGGCTACAAAGTCTTGTGTGGTCTGGGCCAC
ATCCTGTGCCAGGCTCCCTTTTATTGTCTATATCCAGTCAAGCTGGCCTTCAAAATTCAGCACATTG
GGTCTTTCTGCCACAGGCCCTTGACATACATGCTTTTCCCTCTGTACATGTGATTTTCTTTCATCCCC
GCTTTCCAAAGCTAACTCTGACTCATCCTTCAGACTGCATCTTAATGATTACTTCCTTCTGGAAGACTC
CCCTAATTTTCAACTAGATGAGCTCTCCCTTATATGTTGAAAACAAGTATTTCTTGTGTGAACAGAA
ATGTTTAAATGTCTGCCTTCCCACTCAGCCACAAGCTCTGTGAGGATGGGGACTCAGTCTTTCATCCC
TCAGAGCCTAACATAGTCTCAGTTAATACGTGATAGCTGGAGCGTGACATTTTCATTAGGAACAGTGAGG
CTAGTTGAGAAATTATACCTCACCAGTATATAGGAGAGTACCATTAGGAATTCCTAATGGGATGAAAC
CTGAATTCATGTCTAGACCTGCCACCTGCTGGCTGTGCTGTTTTAAATAAATTGCAGATCCTCTCAGA
GCCTGTTTCTCATATATAAAATAATGGGGATCCTAACAGCCCCCTCATATGGTTGTAGCAGGAATTTA
ATTAAATAATATATGAAAAGCATTTTTTAACTGTAACCACATACATAAATATCCTTTGTCACTTTTC
CCAGTGATGACACACTGGGTCACTAAACATCTATTGGATGGGTG
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152749.3</u>
Locus ID:	222255
MW:	45.8