

Product datasheet for SC213039

PELI3 (NM_001098510) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PELI3 (NM_001098510) Human 3' UTR Clone
Symbol:	PELI3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001098510
Insert Size:	1200 bp
Insert Sequence:	>SC213039 3'UTR clone of NM_001098510 The sequence shown below is from the reference sequence of NM_001098510. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CGCCTCATTTTCCAGGGCCCGTGGATTAGGCTCCCTGGGGCCCCCTGCTGCTGTGCCACCTGCCAC
CCAGGTCCCCACCTCCTGCAGCCAGAGGGAGCTCTGCATGTGGGACACTCCCTGCTGGCACAGCCACA
CCAGATGGACATGTTGGATGGGCTGTGCCCTTCCCCCAACTGTGGCCCCCAAGGAGGTCCCCAAGAT
CTCCACCCAGTCATGCTGATGGGGCCTTCAGCACCAGCTCTGTCTGGGTCGATGGAGGAAAGCCAG
CCCCATGGCCTTGCCCTTCTGGGGCATCCACATCGTGCCGCCGACACTGTGTGCCCTGGGGGAGTG
AAGGGGCCAGGGGCTTTGACCCCAAGCTTAGGCTGGCTATGCCCTGAGCCTGCCAACCAGTTTCAGA
CACTCATCCTTCATCCTTCTGCTGCCGCCCTGGGTTCCTCTATAAACCTCGCTGCTCAGCTGCCCCCA
CAAGTCCCACATGCCCGCATGCGTGCAATGCCTCCAGCCCCAAGTGAGTGGTGGGTGGGTGGCAGGC
TGGGCCTAGCACTTGCTGATAGCCACTATGGGCTCGACTCTGTGCTAAGTGGTATCTGGGACACAGAG
GCAATCAGACCTGGTTCCCCCCTGGTGGAGTTCACAGTCTAGTGGAGGACAGACTGTAAACAAATGGC
AACACGGCTTCATGCGTGCCAGGATCAGAGTACAAATCAGGCTGTGGCAGGCCTCACTGGTAATGGA
AGTGAGGGAGGCTTCACGGAGGAGGTGTCATAGGAGCTCAGCCTCAACGGCTGCTCGGAACCCTCCGAG
TCCACAAGATGGGAAGCAAGTCCACAAGATGGGAAGCGTGCTCTATACAGAAGGGAGGGCTGTGCACA
GGCTCAGGTGTCAGTAAGCCAGAAATCCAGGGACAGAGTCAGGATGGCTGGAGCATGGACTGCGGGGGG
AACAGAGAAGGGGGTGGAGCAAAAGGTTGGGACCATGTGACCAGGGAAGATGCTGGATGGAGGATTCTG
AACCTCATTTGTGGACGGTCAGGAGTCAACAAAGGGAGGAGAGGACTCACTGGTCAAGAGTGTGCGC
TGTGACACAAGGGCCATCATCCCTGAATGCCACCATGTGCCCGCACTGTGCTGATGATAATAAT
AACTACGTCCATTGAGCACTCTCCAGA
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001098510.2
Summary:	The protein encoded by this gene is a scaffold protein and an intermediate signaling protein in the innate immune response pathway. The encoded protein helps transmit the immune response signal from Toll-like receptors to IRAK1/TRAF6 complexes. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2011]
Locus ID:	246330
MW:	41.8