

## Product datasheet for **SC213041**

### ZNF3 (NM\_032924) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** ZNF3 (NM\_032924) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** ZNF3  
**Synonyms:** A8-51; HF.12; KOX25; PP838; Zfp113  
**ACCN:** NM\_032924  
**Insert Size:** 1200 bp  
**Insert Sequence:** >SC213041 3'UTR clone of NM\_032924

The sequence shown below is from the reference sequence of NM\_032924. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ACCGAGTTAAATATCAGAGAGTCCACGTGAAAGAGCCACACACCCATTTTCCTCACTTTCCCTGAGTCT
CAAGAGCTCTTGCCCTTACCCTATAAATCTCAACAGCTTAGGATGTGTCCCTTTCAACTCAGACTTTTCA
TTTTAGAGAATGGGGCAGATGGGGCAAATCGTTGAATTTTCCAGAAATCACACCAAGCCTTAGAAAAGCG
TCAAGGCAAGTGGATGGCGTGCTGGGAATAGAAAGCAGCTCTGGGACCAGTTACCCCATTTAGGAAAGG
AGTTTGCCTAACTGTTTTTCTCACAGCAGAGGAACCTTTCCAAGGTGGGGATGGAAGCACAGTCGGG
ACAGAATTCGATGGATTCCTTTAGTTGGAGTCCGCGTAGTCAGCACAGACAGCAGTGAAAGGACGCTC
TGGGTCCTGTTATTTGCTAGGGAGGGTAAAGGGAGACTATTTCAAAGCTACTGTTCTAGTCCAGCTTT
AAGTTTCGGTAAGAAACATGCTGTTTTGTTTCATGATTTCTGTTAATTATGGAAATTTGGCATTGAGGGA
TTATTTTATTGAGGGTAGAAGAGATTCCAGAATCATCATCTGTGATGATGGTGTCTTTAGGGCTCTTG
GAGCAGCCAGACCATGTTTCCAAGAGAAACCTGGTGATATTGCCAGCAGACCCCTGCCATCCCCCCCA
GTTGTCCTGGGGCTGAATGGGCAAATCTGTCCAACAGCTAGTAACCGGCTGTGAGGGAGAGGGTCAAG
AGCACTTAGCGTTGGCCTCTGATTGCTGTCTCTCTTGTCTCTTCCCACTCCAATGATGAAATGATT
TTCTCTAAATGCCTGGTAAGGATGCTTTCAAGGAGCTCACTTGGCCTGCTTGCCCTGCCCTCTCACCT
CTGACACCCAGCCCCAGGAGCCAGACCACTCTGCCTCCACCTCTGACTTTCAGCAGCTGAAGATTAA
TGCAGAGAAAGAGCAAAGCCCAAAGGGGAGAAAAAACAACAAAAACGAAGCCAGAGA
AGGCAGGCTGTGCTGAGCCTGGTAGCAGAGCTCGTTGCTGGGTGAGGATGGCACTCCCCGAGTTTTCC
CAACCGACGGGATGACCTTCCATTGTGTTTTCCAACCTGCTCATCTGTAGCATAAATGAAGTGAC
TATTAACAGAATAGTTTTTCAGAAGA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_032924.5</a>
<b>Summary:</b>	Involved in cell differentiation and/or proliferation.[UniProtKB/Swiss-Prot Function]
<b>Locus ID:</b>	7551
<b>MW:</b>	45