

Product datasheet for SC213042

KIAA0802 (MTCL1) (NM_015210) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KIAA0802 (MTCL1) (NM_015210) Human 3' UTR Clone
Symbol:	KIAA0802
Synonyms:	CCDC165; KIAA0802; SOGA2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015210
Insert Size:	1220 bp
Insert Sequence:	>SC213042 3'UTR clone of NM_015210 The sequence shown below is from the reference sequence of NM_015210. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTCTTGCTAACTGCCCCCTGGGACTCTAGCCCTGCCCGCCTCACGCTGAACCTGTTCTGCACTA
GCTCCATCCCTAGAGCCCTGCTTCTCCAGGCCGAGAGACCAGCAAAACCGTCGCCCTCCGTCCTGTTGG
GCCCCACATTCCCCACTGCCTCACAGCCTCAGTCACCCGGAGACCCGACGTCCTTGGAGGAGCATGGT
GGCGAGGAGCCGCCGAGGAGCAGCCACACCGAGATGCAAGCTTGATGGATTATCACAGTATAATTCA
CTGTAATTTGCATAACCACACCATCACCATGAACAAAACCTCTGCCCAACAGGAGAGATCTAGTTTTCTC
AAGGTCAAAGAATGTTTTTAAAAACACAAAGCTGCTGAATGTTCAACCTGTGAACTGAGATGTTTCT
AGAATGAAACAGTAAATGTGCCTGTAATAACTTAATTTTTTTCATAGCTCAGAAAACCTATTTTTGTCTC
CATCTTTTTTACACACAGTATATTAACGAAAAGGTAAATAAGGTATAAATAGATTTAAAAAATAAAG
TTTTAAAAATGTACATTTTAAGAGATTCTGAACACCCTCGCTGTCAATACCTGACTGCCTCTGTTAAA
TTTGCAGTGTACATTTTGGTTTCAGTTTATTTCCATGTTGAATTAGAGTGGATTAAGTTAATTTTATTT
TGTCAGTGTACTGTTTTTACGAATTTTTTAATGCTTCAGACTGTCTGATTGAGTGAACCTTTTGTAG
TGAAAAAGCCATGAAGCCAGTAGACAAGACAGATATTCTGTATGCTGGAGGGGATACAGGATGATTTTG
AAAAGGTACAAAGTCCTCAGTGGGCTTAGAAAATTCAGTGTATGATCCTTATATTATCCTACTTGGCTT
GCACGTCTTCGGGTGCATGTATATACCGCTACTGTGCTCGCCATCACCTAAATGTGACTCAGTCTGT
TCCACTGTAATATGTTGTGAATTTCTTGTACTGTACTTTTATTGTTGGTCTTCTGATCGATGATCC
AACAGCAACACCATTTTTAAATTATTGTGAAAAGATTAAGTGGCAATGTACAGAGTTTACTCAAAGTTT
TCTTAAGGAAAACACTACAAAAGTCACAAGGATACCAATGGAACACATGATGATGCCTCTGGGTC
TGTATGAGACCGTGTGAGTAGAAATAAAGCCCTTCTGAGATGGCA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015210.4</u>
Summary:	Microtubule-associated factor involved in the late phase of epithelial polarization and microtubule dynamics regulation. Plays a role in the development and maintenance of non-centrosomal microtubule bundles at the lateral membrane in polarized epithelial cells. [UniProtKB/Swiss-Prot Function]
Locus ID:	23255
MW:	47.4