

Product datasheet for SC213051

EPC2 (NM_015630) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	EPC2 (NM_015630) Human 3' UTR Clone
Symbol:	EPC2
Synonyms:	EPC-LIKE
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015630
Insert Size:	1222 bp
Insert Sequence:	>SC213051 3'UTR clone of NM_015630

The sequence shown below is from the reference sequence of NM_015630. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGC
GAGACAACAGTAGCTATGGAAGTGACATAACCTAAACACGTGGCTCTGACCTGTGCTGATGGTGTGCA
GTCATTTCATATCCAGCTGAATGCAAAAGGCAACACTCTGTGGATCACAGAGTGTAAACAATGGACCTAA
ATGGACTATAGTATTTGGATGTAAATCCATATATGATGTATATTTTGTAAATTGGGAAATCACTA
CCTTGTAAGTATTTTGTATCATCAATATTTTCTGTTACTTGAATAGTAGATATTCATCATCA
TGCTTTTGCACCTTGAATTTGCAACTGAATGGATTTTAAAAATAATTCTTTAATGGGATCATGAGCATG
AAATGGGATCCTGCATCACTTGTTTAACTATTTATTTTGGCATGTTTACATTTTGTATCTTGTAAGAA
TAAATCCAACCTTTGTGTCTAAAAAGTTAAAGATTCATAGCTAGGAAATGAAATCTTGTAATTTTTTC
TAAAGGAAGTGAAGTTTTCACTTGGTTCATTTTGTTCACAATTTGACTAGATGGACTTTTTGGTAA
ATACTTTAGTGGCATTTCACGTGCAATATGAAGTTCAAGGCAAAATAGTATTTTCTATTACTGTGCAG
GGGAAAGGGATGGATCGATACATGCAAAATTAATGTAGTAACTCACTTTTCCATATATTTGAATGTAT
ATTTCTATTTATGATACCAATTTATAAAAAATAATTACACAGAAAAATGGAATAGGAAAAATTATGCA
TCTAGCACATTTAACTGTGCAATATGAAAATTTTCGAGGATTACATTTTATCTGAAGGCTGCATAT
TTTAACTGGCTTTAAACTGTAAACACATCACATAAAAGATACTTTACCAGGTATGTATTGCATTATATC
ATTGCAATAATTATTGGAAGTCTAGATATCGAGCCATCCAGGTGTTGGCGGGGGAGGGTTGTGGCA
AGATTGTCTTTTCAATTTTGGAGAGTTTCTGTGGCTACAAGGCAAGTAACGGGTTGAAAAAGTCTG
ACTGTAAGCGTTGGACACCTTCATAGTGTAGTGTTTTAGTGACTTTTTTATACGGTTCTTGTAATTA
GATACGTGTAGTGGTGTTCAGAAATGTTTGTATGCACTAGTTCAGACAACCTTCCCTGTTACTTGT
CTTGATAAGTGAAACTGCAGGGAAATAAAAAATACATATCAAAACATG
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015630.4</u>
Summary:	May play a role in transcription or DNA repair.[UniProtKB/Swiss-Prot Function]
Locus ID:	26122
MW:	48.1