

Product datasheet for **SC213059**

B3GNT9 (NM_033309) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	B3GNT9 (NM_033309) Human 3' UTR Clone
Symbol:	B3GNT9
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_033309
Insert Size:	1203 bp
Insert Sequence:	>SC213059 3'UTR clone of NM_033309 The sequence shown below is from the reference sequence of NM_033309. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GCTGCAGGCCCTTCCAATGGGACTCCTAGCTCCCCACTACAGCCCCAAGCTCCTAACTCAGACCCAGA
ATGGAGCCGGTTTCCCAGATTATTGCCGTGTATGTGGTTCTTCCCTGATCACCAGGTGCCTGTCTCCAC
AGGATCCCAGGGGATGGGGTTAAGCTTGGCTCCTGGCGGTCCACCCTGCTGGAACCAAGTTGAAACCCG
TGTAATGGTGACCCTTTGAGCGAGCCAAGGCTGGGTGGTAGATGACCATCTCTTGCCAACAGGTCCCA
GAGCAGTGGATATGTCTGGTCCTCTAGTAGCACAGAGGTGTGTTCTGGTGTGGTGGCAGGGAATTAGG
GAATCCTACCACTCTGCTGGATTGGAACCCCTAGGCTGACGCGGACGTATGCAGAGGCTCTCAAGGC
CAGGCCCCACAGGGAGGTGGAGGGGCTCCGGCCGCCACAGCCTGAATTCATGAACCTGGCAGGCACTTT
GCCATAGCTCATCTGAAACAGATATTATGCTTCCCACAACCTCTCCTGGGCCAGGTGTGGCTGAGCA
CCAGGGATGGAGCCACACATAAGGGACAAATGAGTGCACGGTCTACCTAGTCTTTCCTCACCTCCTGA
ACTCACACAACAATGCCAGTCTCCCACTGGAGGCTGTATCCCCTCAGAGGAGCCAAGGAATGTCTTCCC
CTGAGATGCCACCACTATTAATTTCCCCTATGCTTCAACCAACCCCTTGCTCAAAAAACCAATACCCA
CACTTACCTTAATACAAACATCCCAGCAACAGCACATGGCAGGCCATTGCTGAGGGCACAGGTGCTTTA
TTGGAGAGGGGATGTGGGCAGGGGATAAGGAAGGTTCCCCCATTCCAGGAGGATGGGAACAGTCCTGGC
TGCCCCCTGACAGTGGGGATATGCAAGGGGCTCTGGCCAGGCCACAGTCCAAATGGGAAGACACAGTCA
GTCACAAAAGTCGGGAGCGCCACACAAACCTGGCTATAAGGCCCAGGAACCATATAGGAGCCTGAGACA
GGTCCCCTGCACATTCATCATTAACTATACAGGATGAGGCTGTACATGAGTTAATTACAAAAGAGTCA
TATTTACAAAATCTGTACACATTTGAAAACTCACAAAATTGTCATCTATGTATCACAAGTTGCTA
GACCAAAATATTAATAATGGGATAAAATTA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTCGATTCCACGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_033309.3</u>
Locus ID:	84752
MW:	43.7