

Product datasheet for **SC213076**

SLC16A10 (NM_018593) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC16A10 (NM_018593) Human 3' UTR Clone
Symbol:	SLC16A10
Synonyms:	MCT10; PRO0813; TAT1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018593
Insert Size:	1206 bp



[View online »](#)

Insert Sequence:	>SC213076 3' UTR clone of NM_018593 The sequence shown below is from the reference sequence of NM_018593. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon CAATTGGCAGAGCTCAGAATTCAAGCGATCGC CATCTGGAATGTTCAAGAAAGAATCTGACTCTATTATTAAATATCTTACATACCTCCACCAGACTGGACT TGCTTTTTGAATTTAAGCAAGTTTCCTTTCTTTTATACAAATTGCAAATTTTCATATTTTTTAAATCAC ATCCTAGGAATAGCACATAATTGGGAAATAGAACCTTATCACTAGAAGAACCATTTTCTGCCACTAA TATCTCTGATGTTCCATGAGTCTGAGGGCAGAGACTCTGGTATATGAAAACGCTGAAAGTCACATATT GTGAAATTTGAAGCTATCTCAGTAAAAAGCAGCTTTGGAACTGTGAATGATCTTTAGCTTGACAAAT GTTTAAAAATACCTCAGGCTATACTGAAAGGGTGCAGTTTGGTATAGGAGTGGAAATATTTGTTTGTTA ATGATGTCTTCAGTTCTGGTACCTCTGTTTTACTTTCTTATGCTCTTTGGAACTTTTGC AAAATTTAA GCCTGGGTCTAGATAATACCAGATCTACCTAAACCTCAAGTCTATGTTAAAGTTGCTTTCTGCTGTTA AATAAGCTATGATATTAAGATATTCTGACTTGCTCCAGTGTCAAGGGACCTTCTGGGAGCAGGTGCTAAC ATAGTGTTCAGAAATCAATATGTGAGATGAAAAGGATCCCTCCAGGAGGATCCTGAGCTGTTAGAAATC ATTTAAGTTTACAGCGTTGTTCCCTTTGCGTTTGCAGTGCCTTTTACTCAAGTAGCCAGAAACACCCAC GTTTCTGAATTTGTTTAACTGTAACAATAAAGTAAATAGATGCATGAAAGATATTCTGGCGATTGTA ACTTAGAATTTTCTGACTTCTGGATTGTTGGCACTAGAACCTGATATTTAAAGTCTTACTGAGCAGCT ATCAAGTGGCAGTTACAGGCACAAATTTGGTGGAGGCTGGAGGATGGGGAGGGGAGCAAAACCTTTATAT TTGTGAAGAAAATATCTGTAGCTGATAGAAATAATTGCTTAAATTGGTTTATGAAATTAATGAGTCTGAA AAGGTTAAAAGCACTTATAAAAAGAACCAAGTCTACATTTCCAGAACTTCTGGCAAAAATTTGCACTC ATATTATTTATCTATGAACATTCCATTGTTTTTTTTTGCTATTTATATACAGATTATCATAAGAAAGC TCTCAGTTTGAGGACC ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_018593.4
Summary:	SLC16A10 is a member of a family of plasma membrane amino acid transporters that mediate the Na(+)-independent transport of aromatic amino acids across the plasma membrane.[supplied by OMIM, Apr 2004]
Locus ID:	117247