

Product datasheet for **SC213084**

CLDND1 (NM_001040200) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CLDND1 (NM_001040200) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CLDND1
Synonyms: C3orf4; GENX-3745; Z38
ACCN: NM_001040200
Insert Size: 1206 bp
Insert Sequence: >SC213084 3'UTR clone of NM_001040200
 The sequence shown below is from the reference sequence of NM_001040200. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
ACCTTAATGAAGGCATATCGTGTGGCATGAGCAAGAAACTGCCTGCTTTACAATTGCCATTTTTATTTT
TTTTAAATAAATACTGATATTTTCCCCACCTCTCAATTGTTTTAATTTTTATTTGTGGATATACCATT
TATTATGAAAATCTATTTTATTTATACACATTCACCACTAAATACACACTTAATACCACTAAAATTTAT
GTGGTTTACTTTAAGCGATGCCATCTTTCAAATAAACTAATCTAGGTCTAGACAGAAAGAAATGGATAG
AGACTTGACACAAATTTATGAAAGAAAATTGGGAGTAGGAATGTGACCGAAAACAAGTTGTGCTAATGT
CTGTTAGACTTTTCAGTAAACTAAAGTAACTGTATCTGTTCACTAAAACTCTATATTAGTTTCTTT
GGGAAACCTCTCATCGTCAAACTTTATGTTCACTTTGCTGTTGTAGATAGCCAGTCAACCAGCAGTAT
TAGTGCTGTTTTCAAAGATTTAAGCTCTATAAAATTTGGGAAATTATCTAAGATCATTTTCCCTAAGCAT
TGACACATAGCTTCATCTGAGGTGAGATATGGCAGCTGTTTGTATCTGCACTGTGTCTGTCTACAAAA
GTGAAAAATACAGTGTTTACTTGAAATTTAACTTTGTAAGTGAAGTCCAGTTCAGCCGGGGCGAG
GATTAGTATTATTTTTAACTCTCCGTAAGATTTTCAGTACCACCAAAATTGTTTGGATTTTTTTCTTT
CCTCTTACATACCAGGTTTATTAAGGTGTGCTTTCTTTTTACATTATATTACAGTTACAAGGTAAAA
TTCTCAACTGCTATTTATTTATCCAGCCAGTACTATAAAGAACGTTTACCATAATGACCTCCAG
AGCTGGGAAACCTACCACAAGATCTAAAGTTCTGGCTGTCCATTAACCTCCAATATGGTCTTTATTTT
TTGTGGTAATATGATGTGCCTTTCTTGCTAAATCCCTTCTGGTGTGTATCAACATTTTAAATGTC
TTCTAATTCAGTCATTTTTTATAAGTATGTCTATAAACATTGAACCTTAAAAAACTTATTTATTTATT
CCACTACTGTAGCAATTGACAGATTAATAAATGTAATTCATAATTTCTTACCATAACCTCAATGTCT
TTTTTAAAAAATAAAATTAATAAATGAAAAGAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001040200.2</u>
Locus ID:	56650
MW:	47.8