

## Product datasheet for **SC213086**

### CLDND1 (NM\_001040181) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** CLDND1 (NM\_001040181) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** CLDND1  
**Synonyms:** C3orf4; GENX-3745; Z38  
**ACCN:** NM\_001040181  
**Insert Size:** 1206 bp  
**Insert Sequence:** >SC213086 3'UTR clone of NM\_001040181  
 The sequence shown below is from the reference sequence of NM\_001040181. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
ACCTTAATGAAGGCATATCGTGTGGCATGAGCAAGAAACTGCCTGCTTTACAATTGCCATTTTTATTTT
TTTTAAATAAATACTGATATTTTCCACCTCTCAATTGTTTTAATTTTTATTTGTGGATATACCATTT
TATTATGAAAATCTATTTTATTTATACACATTCACCACTAAATACACACTTAATACCACTAAAATTTAT
GTGGTTTACTTTAAGCGATGCCATCTTTCAAATAAACTAATCTAGGTCTAGACAGAAAGAAATGGATAG
AGACTTGACACAAATTTATGAAAGAAAATTGGGAGTAGGAATGTGACCGAAAACAAGTTGTGCTAATGT
CTGTTAGACTTTTCAGTAAACTAAAGTAACTGTATCTGTTCACTAAAACTCTATATTAGTTTCTTT
GGGAAACCTCTCATCGTCAAACTTTATGTTCACTTTGCTGTTGTAGATAGCCAGTCAACCAGCAGTAT
TAGTGCTGTTTTCAAAGATTTAAGCTCTATAAAATTTGGGAAATTATCTAAGATCATTTTCCCTAAGCAT
TGACACATAGCTTCATCTGAGGTGAGATATGGCAGCTGTTTGTATCTGCACTGTGTCTGTCTACAAAA
GTGAAAAATACAGTGTTTACTTGAAATTTAACTTTGTAAGTGAAGTCCAGTTCAGCCGGGGGAG
GATTAGTATTATTTTTAACTCTCCGTAAGATTTTCAAGTACCACCAAAATTGTTTGGATTTTTTTCTTT
CCTCTTACATACCAGGTTATTTAAAGTGTGCTTTCTTTTTACATTATATTACAGTTACAAGGTAAAA
TTCTCAACTGCTATTTATTTATCCAGCCAGTACTATAAAGAACGTTTACCATAATGACCCCTCCAG
AGCTGGGAAACCTACCACAAGATCTAAAGTTCTGGCTGTCCATTAACCTCCAACCTATGGTCTTTATTTT
TTGTGGTAATATGATGTGCCTTTCTTGCTAAATCCCTTCTGGTGTGTATCAACATTATTTAATGTC
TTCTAATTCAGTCATTTTTTATAAGTATGTCTATAAACATTGAACCTTAAAAAACTTATTTATTTATT
CCACTACTGTAGCAATTGACAGATTAATAAATGTAACCTCATAATTTCTTACCATAACCTCAATGTCT
TTTTTAAAAAATAAAATTAATAAATGAAAAGAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_001040181.2</a></u>
<b>Locus ID:</b>	56650
<b>MW:</b>	47.8