

Product datasheet for **SC213088**

CLDND1 (NM_001040183) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CLDND1 (NM_001040183) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CLDND1
Synonyms: C3orf4; GENX-3745; Z38
ACCN: NM_001040183
Insert Size: 1206 bp
Insert Sequence: >SC213088 3'UTR clone of NM_001040183
The sequence shown below is from the reference sequence of NM_001040183. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
ACCTTAATGAAGGCATATCGTGTGGCATGAGCAAGAAACTGCCTGCTTTACAATTGCCATTTTTATTTT
TTTTAAATAAATACTGATATTTTCCACCTCTCAATTGTTTTAATTTTTATTTGTGGATATACCATT
TATTATGAAAATCTATTTTATTTATACACATTCACCACTAAATACACACTTAATACCACTAAAATTTAT
GTGGTTTACTTTAAGCGATGCCATCTTTCAAATAAACTAACTAGGTCTAGACAGAAAGAAATGGATAG
AGACTTGACACAAATTTATGAAAGAAAATTGGGAGTAGGAATGTGACCGAAAACAAGTTGTGCTAATGT
CTGTTAGACTTTTCAGTAAAATAAGTAACTGTATCTGTTCACTAAAACTCTATATTAGTTTCTTT
GGGAAACCTCTCATCGTCAAACTTTATGTTCACTTTGCTGTTGTAGATAGCCAGTCAACCAGCAGTAT
TAGTGCTGTTTTCAAAGATTTAAGCTCTATAAAATTTGGGAAATTAAGATCATTTCCTAAGCAT
TGACACATAGCTTCATCTGAGGTGAGATATGGCAGCTGTTTGTATCTGCACTGTGTCTGTCTACAAAA
GTGAAAAATACAGTGTTTACTTGAAATTTAACTTTGTAAGTCAAGAAATCCAGTTCAGCCGGGCGAG
GATTAGTATTATTTTTAACTCTCCGTAAGATTTTCAGTACCACCAAAATGTTTTGGATTTTTTTCTTT
CCTCTTACATACCAGGTTATTTAAAGTGTGCTTTCTTTTTACATTATATTACAGTTACAAGGTAATA
TTCTCAACTGCTATTTATTTATCCAGCCAGTACTATAAAGAACGTTTCACCATAATGACCTCCAG
AGCTGGGAAACCTACCACAAGATCTAAAGTTCTGGCTGTCCATTAACCTCCAATATGGTCTTTATTTT
TTGTGGTAATATGATGTGCCTTTCTTGCTAAATCCCTTCTGGTGTGTATCAACATTTTAAATGTC
TTCTAATTCAGTCATTTTTTATAAGTATGTCTATAAACCATTGACTTTAAAAAATTTATTTATTTATT
CCACTACTGTAGCAATTGACAGATTAATAAATGTAATTCATAATTTCTTACCATAACCTCAATGTCT
TTTTTAAAAAATAAAATTAATAAATGAAAAGAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001040183.2</u>
Locus ID:	56650
MW:	47.8