

Product datasheet for SC213090

DHX40 (NM_024612) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DHX40 (NM_024612) Human 3' UTR Clone
Symbol:	DHX40
Synonyms:	ARG147; DDX40; PAD
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024612
Insert Size:	1203 bp
Insert Sequence:	>SC213090 3'UTR clone of NM_024612 The sequence shown below is from the reference sequence of NM_024612. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CACAGTGACACACGAAAGGAAACAGGCAGGTGGTGAACCTCCAATTCAGGAAGTGGGAAAGGAGC
CAGGAAATGTGCTTCTACTTTGCCAGTTATTTAGACAGCACTACCAAGAGGAGGTGGTCAGCACTTGT
TATTGGCCTATGAACTAAAGCAAATCAAAGCTCATAAATCAAAGCTCATCAGTTCCTATAAATGCAGT
TGTCAAAGAAAAGATTTGGTTGCCATAGTCATAAGCAATGATACATGAAACCAATGAAAGACAGTACAT
GTAATAATATTTTCTCAGTACAATTTTCTGGCCTTAAGTGGTATCAAACGCTGTCAATTGAGATGTTT
TCAAAGAACATTGAGTTGTATTTAATCAGCGTGTACTCCATTTGCATTGAAGCATTAAAAATTTTTTT
CTTAAATCTCTTTAAGGCCCTTCTGTTGCTGTTAGAATAGTGCTATATATCAGGTATGTGACCATTTA
TTTCAGAAGGCTGAACATAAGAGGTTTCTACTCAGCAATACTTAGATGTCTAACTGTTTAAATTGCTACA
GAGCTTTATAGATATTTAGAGAAAAGACTTAATCAATTAGTAAATAAAATTGCCTATGGCAGGATTTCT
TCTTGAATTAATATTAATCCTTAAATTGATTTTTCTGGGATTATACAAATTCCTTTTATATAAAAGTA
TATTGTTTAAACAGTAGCTATAGCCATTAACCAAGGACAGATGATATATATATATATATATATATATATA
TATATATAAGTTCTTTTTAGCTGTACCTACGTACTTATATCAGCACCATGTATGTAGGTGTGATAGTA
CTTTCAAACAGCGCTCCACCTGGCCTACTCTGTTATTTCCACCTGTTTGGGTAGGCGCATTTAACTTC
CATTATGCCAACTTGGGATGGGATTTTGAAGCAGACAACACTATTTTCATCGTGTTTCAAATTGGAAC
CTTGAGGCTAGTTAGTATCACACTCAGGCCACACTCAGCACTTGCCCACTTGTGTTACTGCCTGTAT
TCTAGTTATTTGTGATTTGTCTCCCTCACTAGATTATACGCTCCTTGTGGGCAGGACTGTGCTTTT
TTCATCTTTGTATCTTTTCATGCCTAGCATAGTGCTTTGCACATAGTAGTCACTCAGTGTGTTGTTAAA
TAAAGCTATTAGTGTCAATTAATTAATTAATTAATTAATTAATTAATTAATTAATTAATTAATTAATTAAT
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024612.5</u>
Summary:	This gene encodes a member of the DExH/D box family of ATP-dependent RNA helicases that have an essential role in RNA metabolism. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 17.[provided by RefSeq, Oct 2009]
Locus ID:	79665
MW:	46.3