

Product datasheet for SC213106

CSRP1 (NM_004078) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CSRP1 (NM_004078) Human 3' UTR Clone
Symbol:	CSRP1
Synonyms:	CRP; CRP1; CSRP; CYRP; D1S181E; HEL-141; HEL-S-286
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004078
Insert Size:	1213 bp
Insert Sequence:	>SC213106 3'UTR clone of NM_004078 The sequence shown below is from the reference sequence of NM_004078. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGAGCTGGGGCCTTGGTCCACTCTGAGTAGGCCACCATCACCCACCACACCCTGCCCACTCTGCGCT
TTTCATCGCCATTCCATTCCCAGCAGCTTTGGAGACCTCCAGGATTATTTCTCTGTCAGCCCTGCCACA
TATCACTAATGACTTGAACCTGGGCATCTGGCTCCCTTTGGTTTGGGGTCTGCCTGAGGTCCCACCCC
ACTAAAGGGCTCCCCAGGCTGGGATCTGACACCATCACCAGTAGGAGACCTCAGTGTCTTGGGTCTAG
GTGAGAGCAGGCCCTCTCCCCACACCTCGCCCCACAGAGCTCTGTTCTTAGCCTCTGTCTGCGTGT
CCATCATCAGCTGACCAAGACACCTGAGGACACATCTTGGCACCCAGAGGAGCAGCAGCAACAGGCTGG
AGGGAGAGGGAAGCAAGACCAAGATGAGGAGGGGGAAGGCTGGGTTTTTGGATCTCAGAGATTCTCC
TCTGTGGGAAAGAGGTTGAGCTTCTGCTGTCCTCAGAGTAAGCCTGAGGAGTCCCAGCTTAGGGAGT
CACTATTGGAGGCAGAGAGGCATGCAGGCAGGGTCTAGGAGCCCTGCTTCTCCAGGCCTTGCCTT
TGAGTCTTTGTGAATGGATAGCCTCCCACTAGGACTGGGAGGAGAATAACCCAGGTCTTAAGGACCCC
AAAGTCAGGATGTTGTTTGATCTTCTCAAACATCTAGTTCCTGCTTGATGGGAGGATCCTAATGAAAT
ACCTGAAACATATATTGGCATTATCAATGGCTCAAATCTTCAATTTATCTCTGGCCTTAACCTGGCTC
CTGAGGCTGCGGCCAGCAGAGCCAGGCCAGGGCTCTGTTCTTGCCACACCTGCTTGATCCCTCAGATGT
GGAGGGAGGTAGGCACTGCCTCAGTCTTATCCAAACACCTTTCCCTTTGCCCTGAGACCTCAGAATCT
TCCCTTTAACCACAGACCTGCCTCTTCCACTCCACCCTTCTCCAGGACCCCTAGATCACATCACTCC
ACCCCTGCCAGGCCCCAGGTTAGGAATAGTGGTGGGAGGAAGGGGAAAGGGCTGGGCCTACCGCTCCC
AGCAACTGAAAGGACAACACTATCTGGAGCCACCCACTGAAAGGGCTGCAGGCATGGGCTGTACCAAG
CTGATTTCTCATCTGGTCAATAAAGCTGTTTAGACCAGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_004078.3</u>
Summary:	This gene encodes a member of the cysteine-rich protein (CSRP) family. This gene family includes a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. The LIM/double zinc-finger motif found in this gene product occurs in proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2010]
Locus ID:	1465
MW:	44.3