

Product datasheet for **SC213111**

GPR105 (P2RY14) (NM_001081455) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: GPR105 (P2RY14) (NM_001081455) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: P2RY14
Synonyms: BPR105; GPR105; P2Y14
ACCN: NM_001081455
Insert Size: 1213 bp
Insert Sequence: >SC213111 3'UTR clone of NM_001081455
The sequence shown below is from the reference sequence of NM_001081455. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACAACACTTGAAAGCACAGATACTTTGTAGTTTCCTACCCTCTTCCAAGAAAGACCACGTGTGCATGT
TGTCATCTTCAATTACATAACAGAAATCAATAAGATATGTGCCCTCATCAATAATCATCTCTAGCAC
TGCCATCCAATTTAGTTCAATAAAATTCAAATATAAGTTTCCATGCTTTTTTGTAAACATCAAAGAAAC
ATACCCATCAGTAATTTCTAATACTGACCTTTCTATTCTCTATTAATAAAAAATTAATACATACAAT
TATTCAATTCTATTATATTAATAAAGTTAAAGTTTATAACCACTAGTCTGGTCAGTTAATGTAGAAAT
TTAAATAGTAAATAAACACAACATAATCAAAGACAACCTCACTCAGGCATCTTCTTTCTCTAAATACCA
GAATCTAGTATGTAATTGTTTTCAACACTGTCCTTAAAGACTAAGTTGAAAGCAGGCACAGTTTGATGA
AGGGCTAGAGAGCTGTTGCAATAAAAAGTCAGGTTTTTTCCTGATTTGAAGAAGCAGGAAAAGCTGA
CACCCAGACAATCACTTAAGAAACCCCTTATTGATGATTTTCATGGCACTGCAAAGGAAGAGGAATATT
AATTGTACTTAGCAAGAAAATTTTTTTTTCTGATAGCACTTTGAGGATATTAGATACATGCTAAAT
ATGTTTTCTACAAGACTTACGTATTTAATGAGCCTGGGGTTCTGGTGTAGAAATTTTTTAAGTAGG
CTTACTGAGAGAACTAAATATTGGCATACTGTTATCAGCAACTTCCCCTGTTCAATAGTATGGGAAAA
ATAAGTAGTGGGAAAAAGACACACCCACACCGTAGAACATATATTAATCTACTGGCGAATGGGAAAG
GAGACATTTTCTTAGAAAGCAAATAAAGTTGATTTTTTAAATCTAAAAATTTACATTAATGAGTGCAA
AATAACACATAAAATGAAATTCACACATCACATTTTTCTGGAAAACAGACGGATTTTACTTCTGGAGA
CATGGCATAACGGTTACTGACTTATGAGCTACCAAACTAAATCTTTCTCTGCTATTAAGTGGCTAGAA
GACATTCATCTATTTTTCAAATGTTCTTTCAAACATTTTTATAAGTAATGTTTGTATCTATTTTCATGC
TTACTGTCTATATACTAATAAAGAAATGTTTAAACTG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001081455.2
Summary:	The product of this gene belongs to the family of G-protein coupled receptors, which contains several receptor subtypes with different pharmacological selectivity for various adenosine and uridine nucleotides. This receptor is a P2Y purinergic receptor for UDP-glucose and other UDP-sugars coupled to G-proteins. It has been implicated in extending the known immune system functions of P2Y receptors by participating in the regulation of the stem cell compartment, and it may also play a role in neuroimmune function. Two transcript variants encoding the same protein have been identified for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	9934
MW:	47.9