

Product datasheet for **SC213151**

ZNF407 (NM_017757) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF407 (NM_017757) Human 3' UTR Clone
Symbol:	ZNF407
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017757
Insert Size:	1233 bp
Insert Sequence:	>SC213151 3'UTR clone of NM_017757 The sequence shown below is from the reference sequence of NM_017757. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGAGAAAGCAGCGAACTCCAGGAAGCATGAGACGCGCGGCACCTTTACTCAGCACAGGGCAGGTGTGGG
AAGGTCCAGCTTCGGTGGGGGACCGTGTTCCTGAGCTTCATCTGAAACCTTCAAACCATGAGGACAA
GGCTCCCGTGAGCTCTGAGCATGCCCTCCAGCGAGAGTCACACTGGCCACCAGCCAGGCGCCACAGA
GGGTACCGTGGGCTGGGCTCGGGGAGCAGGCTGCCAAGTGCAGGGGAGGGCGGCGCAGGCCGCACA
GGGAGCTCCGGTCCACTGGGGGCCCTTCGATCAGTGGCCTCCCGCTTCGTCTGGCCGCTGTGCTGAAG
AGAAGCCAAGTGTTGTTGGTGTTTTCTCTCCAAGTGTTTTCCATTTAGTTATCAGAAGGTCATGG
CCGTGGGGAAGTGTTGAAGATACCCCTCCTGGCTTGGGGTGACCTGCTTGTGCAGTCAGCATGTAGC
TGCCTTTCCATTTCACTCTACTGGGCTAAAAATTGCAGCTACAAGTGTACCATCTTGAAGCAGTCC
ACTTCCATTCAATTTTTTTTTTTAATTTTGAATAACAGTGCCCCATACCAAGGAAGCCTGCTAGC
TCATTTTCATGTATAAATTTCCCATCTTCAAACAGTTTAGGTGTATTTGTTGCTCTGGTCACATTCTGCA
TAAAGAAATCCTCTTAAGCCTATGGTTAAGAAAAGCCTTGAAGTTTATATTCAGTTAAATATATGTC
GGTGGAGATAGCCAGTGCTTCTAATTTTGAAGTTAGTTTTCATACAGTAAAGCCTAAATGTGAAACGCACA
CGCTGGAAGATATTGTTCTATCAATATTTGCTTTTTATAACAAGGGTTTGTTTCATATTGATGCCATT
TTTGCAGGATTTCTCGTGATTTCTGTCCATATGAAAATGCTGACATTAACATTAACACATGGAGACC
GTGCCCTGTGGCCCTGCCGTGGTGCCAGCATGGTCTGTGTTTCTTGTGGATTACCTGTGGCCCTGC
TGTGGCCACCAGCATGGTCTGTGTCCTCGTGGATTCACTGCAGCTGTGGATGCGAGTTTCTGTCATAA
TCATTTGTTTCTGATACAATTGTTCTTATTCTTTTCCAAACTGTAATAATCTCCTCCCTCAAATG
CAAAGGTTGTTTTGTTCTGTTTCTGTTTCTTTGAAATAAAATTATAACGTTAAAGAT
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_017757.3
Summary:	This gene encodes a zinc finger protein whose exact function is not known. It may be involved in transcriptional regulation. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]
Locus ID:	55628
MW:	45.9