

Product datasheet for **SC213157**

GRAMD2B (NM_023927) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GRAMD2B (NM_023927) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	GRAMD2B
Synonyms:	GRAMD3; NS3TP2
ACCN:	NM_023927
Insert Size:	1438 bp



[View online »](#)

Insert Sequence: >SC213157 3'UTR clone of NM_023927
The sequence shown below is from the reference sequence of NM_023927. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTACAGAAGTTGCTTGAGAATGGTGACTGATCGACCAGATTGCTTGGCCATCGGAATACCTCATGTTT
CCCTTTGAAGAAGGTGCTTCTGAGGCGTTTTGTTTGGAGTGCACCCTGCTGGTCAGAGGTGCAAGCAGA
TGAGAATCCAGACATTGCATGTGGAGGTCTCCAGCTCAGGAAAGTGAGGAGGGAAATAATTTTGGTTCT
TGTGCAATAAAAAGTTGACCTTGACTCTCTGAGGAAGATTTTGGCTTTTTGCCTGAAGAAAAACAGACCC
ATCTCTGGAGGTCTCAGGAAGGGCCAGCGAACACTCTCTGGATAATTACCACGATGGCGTCAGCA
AACACTCCACCCTGTGCCTTTTTAGTCCTTCCCGCCCTCTGCCTCTCCCTTACACCCTCTTAACGAC
TTTCAAATAAAGGATACATCATATACTGACAACTCAATGTGGTCTTTCAAGAATTAGCCATGAGTC
TCAAAAAGGCAATAAATGGCTCTAAGTGGACAGGTTTGGCTTCAAACAAGTAACATCTACATTTTGTCTT
TTTTTTTTTCAGTTCTCCTGTTATGTTCTGGTTGAAATCACCTGTGTGCTTAATTTCTCAATTCCTTTT
TGGCAAGAATATCAAGCAAGGTGAATTTAACATTATGTTTATGTTTGTGTTTGTGCTGTAACATAATAG
TTAATTGGACTGATTCTTACCCAGCCCTGGTCAAGAATCTGTGAGGCATGTGACTGAAGTACTAAATTA
AACTTATTTTGAACCAACCTAATTTTTAAGCCAAAAGGTGTAATAGTGATTTAATACAGGATGAAAA
ACACTGAATTTTTAAGACTGTAGGTGGACTATGTTAGTAGTTTTCAAGCAGGATGTCTGTATTACAGCAT
TCAATAATGCTAAAATCCCTTTCAGCATGAAATTTGTATGTTTTATCCTTTGCTGACTAAAATAAAT
AACTGGTGGTTTTGCTAACTATTTACCATATGGTTGGCTGGTTTTTATTAATAATTGTTTAGGGTATCA
TAAGATCTGTAAGTACAAAGGATATTCTGTTTCTATCTGCAATGCATTTGATAATCATTCTATGAAAT
GTATTTTATTTAATCAGATAAGCTAATAAGCGAATTGGTTACATCGTTTGTATCTGTTGGCCTAGTGGA
CAGCTGTGCCTGGTGCACCTCTGGCTTTTAATAAATACTATTGGTTGAAACACCCTGCTAACACTGTG
TTTTGTGAGTTTCGTAATTCCTTTAACCAAGACTTAACATGCTTTTTAACTCTTTTTCCACCTGCCAT
CTACTACTACCAACTCACTTATACTTCTTCTCCCAAGTAATCATTTCAAACCTAAAAAAGGATTTAA
GTTCTTGGGTTACAACGGTTTAAAGGAAAGAAAATGAACTTTGAAGTCAAAAAAAAAA
ACGCGTAAAGCGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_023927.4](#)

Locus ID: 65983

MW: 55.4